





# ELLIPSE TECHNOLOGIES

## MAGEC™ TRAINING PROGRAM



# The Company



- Founded August 2005
- Focused on the Development of Non-Invasive, Remote Controlled Orthopedic Devices
- 30 Full Time Employees
- Sales coverage in 26 countries worldwide
- CE Mark registration
- 2,500 sq. ft. Facility with Lab Located in Irvine, CA

# MAGEC™ Focus: Early Onset Scoliosis (EOS)




- EOS is a form of Idiopathic Scoliosis
  - Cause is "unknown"
  - Occurs in children before 5yrs in age, predominantly in boys
- EOS can progress rapidly into severe, complex curves that impact thoracic organs, lung capacity and respiratory functions
  - Can be fatal
- Can be treated non-invasively or surgically, depending upon the prognosis



# Anatomy of the spine



- The normal thoracolumbar spine is relatively straight in the coronal plane and has a double curve in the sagittal plane.
- The thoracic spine is convex posteriorly (kyphosis) and the lumbar spine is convex anteriorly (lordosis).
- Normally there should be no lateral curvature of the spine.



# Scoliosis



- Scoliosis is a complicated deformity that is characterized by both lateral curvature and vertebral rotation
- As the disease progresses, the vertebrae and spinous processes in the area of the major curve rotate toward the concavity of the curve.
  - On the concave side of the curve, the ribs are close together. On the convex side, they are widely separated.



# Scoliosis



- As the vertebral bodies rotate, the spinous processes deviate more and more to the concave side and the ribs follow the rotation of the vertebrae.
  - The posterior ribs on the convex side are pushed posteriorly, causing the characteristic rib hump seen in thoracic scoliosis.
  - The anterior ribs on the concave side are pushed anteriorly.

# MAGEC™ Product Description



The Ellipse Technologies, Inc. (ETI) **MAGEC™ Spinal Bracing and Distraction System** is comprised of one or two sterile spinal rods that are surgically implanted using commercially available fixation components (i.e. pedicle screws, hooks and/or connectors).

1. *The system includes a non-sterile hand held External Remote Controller that is used at various times after implant to non-invasively lengthen or shorten the implanted spinal rods.*

2. *The rods include a small internal magnet which allows the rods to be lengthened by use of the External Remote Controller.*



## MAGEC™ Indications for Use



- The implanted rod is used to brace the spine during growth to minimize the progression of scoliosis.



## MAGEC™ Contraindications



- Infection or Pathologic conditions of bone such as osteopenia which would impair the ability to securely fix the device.
- Metal allergies and sensitivities.
- Patient with pacemaker.
- Patient requiring MRI imaging during the expected period device will be implanted

## MAGEC™ Contraindications



- Patients younger than 2 years old.
- Patients weighting less than 25 lb. (11.4 kg)
- Patients and/or families unwilling or incapable of following postoperative care instructions.

## MAGEC™ Warnings



- The ETI MAGEC Bracing and Spinal System Implants are supplied sterile and are for single use only and cannot be reused or resterilized.
- Do not use if the sterile pouch has been damaged or is open.
- Metallic implants can loosen, fracture, corrode, migrate, or cause pain.

## MAGEC™ Selected Precautions



- Do not use this device without proper training in both device implantation and adjustment.
- Assure that distraction length is assessed by X-ray imaging immediately after non-invasive adjustment procedure, and also at a minimum of once every six months.
- Assure that patient with implanted device does not enter MRI unit.

## MAGEC™ Selected Precautions



- During period of implant, if brace is used on patient, brace should not have any magnetic metallic components (steel, etc.) which may affect the implanted magnet.
- During period of implant, patient should not participate in contact or severe sports such as weightlifting, tumbling, gymnastics, rowing, or other high risk activities.
- During period of implant, patient should limit backpack weight to 20 lb. (9 kg) or less.

# MAGEC™ SYSTEM



- Titanium Implant




- External Remote Controller (ERC)



**Rotational Force** 

**Axial Force**

Rare Earth Magnets  
Neodymium Iron Boron (NdFeB)



The interaction between rare earth magnets  
modifies the length of the MAGEC™  
implant.


# MAGEC™ IMPLANTS




- Available in 4.5mm and 5.5mm sizes.
- Designed to interface between rod and commercially available fixation components, i.e., pedicle screws, hooks, and/or connectors.



# MAGEC™ ACTUATOR



# MAGEC™ rod dimensions




- Proximal rod length 26cm
- Distal rod length 12cm
- Actuator length 9cm
- Actuator diameter is 9mm

Actuator must never be bent




**DO NOT BEND**

# MAGEC™ DUAL RODS




- The longer portion of rod should always point towards the head.

# MAGEC™ DUAL RODS




- When using dual rods in a patient, the actuators should be placed at the same height as each other in relation to caudal and cephalad

# MAGEC™ SINGLE ROD “SL”



- Diameter of actuator SL 10.5mm

- STERILE WRAPPED IMPLANT



# MAGEC™ IMPLANT



- Minimally invasive implant & explant
- 4.5mm and 5.5mm sizes available
- 48mm of lengthening available
- Eliminates multiple surgeries
- Bi-directional axial adjustments
- Single or dual rod constructs
- Meets International Magnetic Safety



- EXTERNAL REMOTE CONTROLLER (ERC2)





ERC2 Proper Patient Positioning



# MAGEC™ ERC DISPLAY



Control Panel

# MAGEC™ ERC CASE



- ✓ Patient returns to Dr. office
- ✓ Non-invasive procedure
- ✓ Minutes vs. hours or days
- ✓ Less mental and physical trauma
- ✓ Fewer infections from multiple surgeries
- ✓ LCD readout
- ✓ ERC is portable

# COMPETITION



- VEPTR – Synthes
- Traditional Growing Rods
- Shilla-Medtronic
- Phenix-France

## Summary



- ✓ Minimally invasive procedure
- ✓ Single or dual rod construct
- ✓ Non-invasive adjustments
- ✓ Less traumatic for patients & family
- ✓ Eliminates need for multiple surgeries
- ✓ Reduces risk of infections
- ✓ Less pain & suffering