iFACTOR™
Peptide Enhanced Bone Graft

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Multi-Center Clinical Evaluation of Combination Anorganic Bovine-Derived Hydroxyapatite Matrix (ABM)/Cell Binding Peptide (P-15) as a Bone Replacement Graft Material in Human Periodontal Osseous Defects: 6-Month Results

Raymond A. Yukawa,* Donald P. Callan,[†] Jack T. Krauser,[‡] Gerald H. Evans,[§] Mary Elizabeth Atchellman-Reidy,[§] Kent Moore,[§] Renee Cruz,[§] and J. Brent Scott[¶]

A SYNTHETIC CELL-BINDING PEPTIDE (P-15) combined with anorganic bovine-derived hydroxyapatite bone matrix (ABM) was compared to demineralized freeze-dried bone allograft (DFDBA) and open flap debridement (DEBR) in human periodontal osseous defects in a controlled, monitored, multi-center trial. Following appropriate initial preparation procedures, flap surgery with defect and root debridement was performed. Three osseous defects per patient were treated randomly with one of three procedures after surgical preparation. Appropriate periodontal maintenance schedules were followed, and at 6 to 7 months re-entry flap surgery was performed for documentation and finalization of treatment. Analysis of variation (ANOVA) and *t* test analyses of patient mean values from 31 patients revealed that the combination ABM/P-15 grafts demonstrated significantly better mean defect fill of 2.8 ± 1.2 mm (72.3%) versus a mean defect fill of 2.0 ± 1.4 mm (51.4%) for defects treated with DFDBA ($P < 0.05$) and a mean defect fill of 1.5 ± 1.3 mm (40.3%) ($P < 0.05$) for defects treated with DEBR. Other hard tissue findings showed similar clinically superior results with the use of ABM/P-15. Relative defect fill results showed 87% positive (50% to 100% defect fill) responses with ABM/P-15, 58% positive responses with DFDBA, and 41% positive responses with DEBR. There were 8 to 9 times more failures (minimal response) with DFDBA and DEBR (20% to 29% frequency) than with ABM/P-15. Soft tissue findings showed no significant differences among treatments except for greater clinical attachment level gain with ABM/P-15 compared to DEBR. These results suggest that the use of the P-15 synthetic cell binding peptide combined with ABM yields better clinical results than either DFDBA or DEBR. Further studies are needed to determine the relative roles of the ABM and/or the P-15 in these improved results. *J Periodontol* 1998;69:655-663.

Key Words: Hydroxyapatite/therapeutic use; bone, freeze-dried; bone, demineralization; cell attachment protein; multi-center studies; bone matrix; comparison study; surgical flaps; periodontal diseases/therapy.

Periodontists have used various bone replacement graft (BRG) materials in attempting to regenerate bone lost as

a result of periodontal disease, but predictable high rates of success have not been achieved. Autogenous bone grafts; allografts such as non-decalcified (FDBA) and decalcified (DFDBA) freeze-dried bone; bio-compatible ceramics such as beta-tricalcium phosphate and synthetic hydroxyapatite; bioactive glasses, polymers; and natural eluted in human intrabony defects generally positive clinical results.¹⁻¹⁹

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*Department of Periodontics, Louisiana State University School of Dentistry, New Orleans, LA.
[†]Private practice, Little Rock, AR.
[‡]Retired, Memphis, TN.
[§]Retired, Memphis, TN.
[¶]Retired, Memphis, TN.

CERAPEDICS
Enhancing the Science of Bone Repair

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Ralph Mobbs

Prince of Wales Hospital
University of New South Wales

INTRODUCTION

Anterior Lumbar Interbody Fusion (ALIF) is a technique utilized for a variety of indications & pathologies, including the management of Degenerative Disc Disease, Spondylolisthesis, Adjacent segment degeneration, Foraminal stenosis and Symptomatic non-union of a posterior fusion.

There are a variety of graft options available to the surgeon including Autograft, Bone Morphogenic Proteins, Allografts and a variety of synthetic based bone graft substitutes (1). This study aims to compare early fusion rates and complication rates of a variety of available graft options, including Autograft, INFUSE and iFactor.

METHODS

Over an 18 month time period, data from 35 consecutive ALIFs has been prospectively collected; including 10 autograft (Iliac crest), 10 Infuse (rhBMP-2) and 15 i-Factor (P-15) grafts. Patients were reviewed at 6, 12, 24 and 52 weeks postoperatively. (NOTE: Clinical outcomes are the subject of a separate study) Early radiological outcomes based on fine cut coronal CT assessment, and graft associated complications were prospectively recorded.

10X2000 bone graft (Autograft) in, Washington, DC 10X2000 bone graft (Autograft) in, Washington, DC 10X2000 bone graft (Autograft) in, Washington, DC

In this study, the bone graft was placed in the anterior lumbar interbody space of the lumbar vertebrae and the bone graft was secured in place by the body's natural healing process.

RESULTS

Early fusion rates were positive throughout all patient groups. By 3 months, there was early bridging bone on coronal CT in 3/10 patients in the autograft group, 7/12 in the Infuse group and 13/15 in the iFactor group. By 12 months there was one non-union in the autograft group (smoker / diabetic) and no non-union in the Infuse group. The average follow-up in the i-Factor group is 5 months, however fusion rates at time of reporting were 14/15 in the i-Factor group. Complication rates at the 3 month mark include; 3/10 patients suffering from significant graft site pain in the autograft group, 1/10 patients suffering a severe BMP related radiculitis, and no complications in the i-Factor group.

CONCLUSIONS

From this early data analysis with regards to fusion and complications, all groups demonstrated early fusion. Bone bridging appears to develop faster in the Infuse and iFactor groups, based on the 3 month data, compared with Autograft. The major graft related complications include residual graft site harvest pain, followed by a severe Infuse radiculitis in one patient. Issues of cost and graft related complications will require longer term follow-up and larger patient cohorts to define more accurately. i-Factor demonstrates very rapid bone bridging/integration with a low complication profile.

i-Factor Results

REFERENCES

1. Chey AA, Mobbs RJ. Bone graft substitutes in anterior cervical discectomy and fusion. *Exp Surg J*. April 18 (448-54) 2000
2. Follmar G, Bardenheir M, Weyand M, Zuo X, Bieger C. Differences in early osteointegration and bone interdigitation in anterior lumbar interbody fusion with INFUSE-2, apogee bone paste, autograft, and autograft. *Bone*. Aug 6(2):287-78. 2000

CONTACT & DISCLOSURE

Dr Ralph Mobbs MD FRACS
Prince of Wales Private Hospital
Sydney Spine Clinic

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© 2000 Mobbs