

The  in Bone Marrow Stimulation

**Joint**  
rep<sup>®</sup>


**OLIGO**  
MEDIC<sup>+</sup>

**Joint+rep<sup>®</sup>**

# Target patients

## The status quo for cartilage repair

- Microfracturing represents the most common cartilage repair technique.<sup>1</sup>
- 63% of arthroscopies reveal the presence of asymptomatic cartilage lesions.<sup>2</sup>



ICRS Grading Scale\*

### Grade I

Cartilage with softening and swelling

### Grade II

Partial thickness defect with fissures on the surface that do not reach subchondral bone or exceed 1.5cm in diameter

### Grade III

Fissuring to the level of subchondral bone in an area with a diameter more than 1.5cm

### Grade IV

Exposed subchondral bone

Bone marrow stimulation  
(BMS) + JointRep<sup>®</sup>

## PROS OF MFX

- Easy to apply with no specific infrastructure needed
- Single-step, easy to apply arthroscopically/mini-open
- Safe with relatively lower cost (when compared to culture cell-based and grafting solutions)

## Could we improve the results of MFX alone?

- Risk of retraction and/or disruption of the blood in the cartilage lesion.<sup>3</sup>
- Inconsistent quality of the resulting healing tissue.<sup>4</sup>
- Eventual deterioration of knee function on the medium/long term.<sup>5</sup>








# Chitosan

## A polyglucosamine polymer

- has strong uses for biomedical applications such as wound healing and is one of the most abundant polysaccharides in nature.
- is a particularly suitable polymer base for a cartilage repair scaffold due to its good biocompatibility and chemical similarity to native joint macromolecules.
- enjoys a unique status as a natural cationic polysaccharide, which confers its bioadhesive property.




The JointRep<sup>®</sup> matrix is composed of completely deacetylated chitosan, also known as polyglucosamine. At physiological pH, chitosan is partially cationic, which results in its bioadhesive property.

“Simplicity is a core guiding principle and a key component of our approach to managing the risk-benefit challenge inherent to developing innovative treatments.”

# Joint+rep<sup>®</sup>


## The + in Bone Marrow Stimulation



Transversal view: Joint+rep<sup>®</sup> protects the blood clot and acts as a bioscaffold for the stromal cells.

### WHAT

- A second generation chitosan-based, easy to use arthroscopically-delivered bioscaffold.<sup>6</sup>
- Designed to be used in conjunction with a single-step bone marrow stimulation procedure.
- A cost-effective first-line treatment for cartilage defect repair.







## WHY

**Simplicity:** Prepared in under one minute, JointRep® is an off-the-shelf solution injected to fill any shape defect (single or multiple lesions) with no wait time after implantation to begin closing.

**Safety:** Isotonic with neutral pH, non-toxic and highly biocompatible, JointRep® has been used in more than 7000 clinical cases as of 2021.

**Efficacy:** Improved consistency in clinical outcomes (vs MFX) through protection of the blood clot.

**Joint**  
**+rep®**


The + in Bone Marrow Stimulation

## HOW

- Once mixed, JointRep<sup>®</sup> forms a flowable viscous gel and can be applied arthroscopically until the defect is filled up to every edge, easily adapting to any lesion shape, size and contour.<sup>7</sup>

In adapting itself to the lesion (and not vice versa), JointRep<sup>®</sup> can preserve healthy cartilage from unneeded debridement.

- Adheres firmly to the lesion surface and surrounding cartilage rim due to its cationic nature.
- Helps to stabilize and protect the BMS generated clot, forming a "hybrid clot" and allowing for optimal healing outcomes in a BMS procedure.



## Clinical Data

Microfractures and hydrogel scaffolds in the treatment of osteochondral knee defects: A clinical and histological evaluation, Pipino et al (2019)

A controlled post-market clinical study (Bologna University, Italy / LUdeS HEI, Switzerland / Stanford University, USA) comparing microfracture with JointRep (test group) versus standard arthroscopic microfracture alone (control group).

A total of 69 patients (46 test, 23 control) aged 18-75 having Outerbridge grade III or IV knee joint cartilage damage were enrolled.

### Patient Inclusion

- No restriction on lesion size
- Associated conditions such as previous partial meniscectomy, cruciate ligament lesions or previous failed microfracture
- 18-75 year old patients

### Primary Endpoint

- Clinical improvement assessed by the Western Ontario and McMaster Universities Arthritis Index (WOMAC) at 6, 12 and 24 months post-op
- Safety assessed by the recording of Adverse Events

### Microfracture Technique


- Standard Steadman arthroscopic microfracture
- No autologous materials mixed with the device
- No special instruments utilised for administration of JointRep

### Rehabilitation Program

- Following arthroscopy, all patients were first allowed to weight bear as tolerated (WBAT) immediately after the surgery; the use of a contralateral cane for 5-7 days postoperatively was suggested
- On day 15 post-op, the patients were allowed to start formal standard physical therapy, including quadriceps electro stimulation, swimming, and the use of a stationary bike for a reduced period of 3 weeks only.



**Improved clinical outcomes for the JointRep Group  
at 6, 12, 24 and 36 months post-operation\***



**\*up to 90%+ improvement in clinical  
outcomes sustained up to 36 months<sup>6,8</sup>**

**Patient Statistics and Defect Size**

Patient Statistics				
	Study Group (JointRep + Microfracture)		Control Group (Microfracture alone)	
	N=46		N=23	
<b>Total number of patients</b>	N=46		N=23	
<b>Patients, age</b>	54.5 ± 9.5 (26 - 72)		56.6 ± 7.6 (44 - 70)	
<b>Patients, gender</b>	Male 29 (63%)	Female 17 (37%)	Male 11 (47.8%)	Female 12 (52.2%)
<b>Treated knee</b>	Right 25 (54.3%)	Left 21 (45.7%)	Right 19 (69.6%)	Left 4 (30.4%)
<b>Grade (outerbridge)</b>	III 10 (23%)	IV 36 (78%)	III 6 (39%)	IV 17 (61%)
<b>Associated lesions</b>	Lesion of Meniscus 98%	Patello-femoral 2%	Leison of Meniscus 100%	Patello-femoral 0%
<b>Previous microfracture</b>	1	0	0	0
<b>Average defect size</b>	2.8 cm <sup>2</sup>	2.5 cm <sup>2</sup>	2.7 cm <sup>2</sup>	2.5 cm <sup>2</sup>



## Case Study

### Patient

- 23 y.o. male
- Isolated full thickness chondral lesion on lateral femoral condyle
- 5 cm<sup>2</sup> after debridement
- JointRep implanted after BMS


### Surgeon

- Dr Mark Clatworthy  
Auckland Bone and Joint Surgery


### 9 Month post-operative MRI

- Complete filling of lesion
- No bone marrow reaction beneath treated lesion

Pre op MRIs



Post op MRIs





## JointRep<sup>®</sup>

- + An off-the-shelf, 2<sup>nd</sup> generation chitosan-based hydrogel prepared in under one minute to complement a single-step bone marrow stimulation procedure that can be injected arthroscopically, filling even odd shaped defects (single or multiple lesions) with no wait time before closing.
- + Helps to stabilize and protect the bone marrow stimulation-generated clot, forming a “hybrid clot” and allowing for optimized tissue healing within a properly prepared articular cartilage defect.
- + Published clinical data demonstrating improved pain and functional outcomes VS microfracture-alone.





# References

1. Medina J, Garcia-Mansilla I, Fabricant P, Kremen T, Sherman S, Jones K. Microfracture for the Treatment of Symptomatic Cartilage Lesions of the Knee: A Survey of International Cartilage Regeneration & Joint Preservation Society. *Cartilage*. 2020;00(0):1-8. doi:10.1177/1947603520954503
2. Widuchowski W, Widuchowski J, Trzaska T. Articular cartilage defects: Study of 25,124 knee arthroscopies. *Knee*. 2007;14(3):177-182. doi:10.1016/j.knee.2007.02.001
3. Furukawa T, Eyre D, Koide S, Glimcher M. Biochemical studies on repair cartilage resurfacing experimental defects in the rabbit knee. *The Journal of Bone & Joint Surgery*. 1980;62(1):79-89. doi:10.2106/00004623-198062010-00012
4. Frisbie D, Oxford J, Southwood L et al. Early Events in Cartilage Repair After Subchondral Bone Microfracture. *Clin Orthop Relat Res*. 2003;407:215-227. doi:10.1097/00003086-200302000-00031
5. Mithoefer K, Williams R, Warren R et al. The Microfracture Technique for the Treatment of Articular Cartilage Lesions in the Knee. *The Journal of Bone & Joint Surgery*. 2005;87(9):1911-1920. doi:10.2106/jbjs.d.02846
6. Pipino G, Risitano S, Alviano F et al. Microfractures and hydrogel scaffolds in the treatment of osteochondral knee defects: A clinical and histological evaluation. *J Clin Orthop Trauma*. 2019;10(1):67-75. doi:10.1016/j.jcot.2018.03.001
7. Please refer to the Instructions for Use (IFU) for a complete list of indications and contraindications
8. Pipino G, Vaccarisi D, Mardones R, Indelli PF. Articular Cartilage Lesion of the Knee: 3 Years Follow-up. 20th EFORT Annual Congress, Lisbon, Portugal. June 6th, 2019. Poster #2378.

**OLIGO MEDIC+** was founded in 2010 by pioneers of chitosan as a material for biomedical applications, with over 25 years of experience developing medical devices

For more visit: [www.oligomedic.com](http://www.oligomedic.com)

OligoMedic Inc  
500 blvd Cartier O  
Laval, QC  
H7V 5B7  
(450) 680-3371

