MEDPRIN

NeoDura

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
Altmetric

Original Research Paper

A novel biomimetic composite substitute of PLLA/gelatin nanofiber membrane for dura repairing

Kunxue Deng , Yaya Yang, Yiquan Ke, Chengyi Luo, Man Liu, Yuting Deng, ...show all

Pages 819-829 | Received 18 Mar 2017, Accepted 24 Jun 2017, Published online: 12 Jul 2017

 Download citation <https://doi.org/10.1080/01616412.2017.1348680> Check for updates Full Article Figures & data References Citations Metrics Reprints & Permissions Get access Выбрать язык

Translator disclaimer

Abstract

Objective: Biomimetic design will significantly improve growth and regeneration of dural cells and tissue for better repairing effects and fewer complications in repairing the native dura. This study designed a novel composite, biomimetic substitute based on the characteristics of native dura extracellular matrix.


Methods and results: This substitute is expected to rapidly induce cell adhesion, migration, and fast regeneration of neotissue. The material characteristics (contact angle, surface charge, and zeta potential were evaluated), in vitro biological characteristics (cell stretch, connections between cells, cell proliferation) and in vivo tissue regeneration capability of this substitute were evaluated, compared to those of collagen dura substitute, the mostly used dura substitute. The results showed that the surface properties of this composite substitute were more biomimetic to native extracellular matrix than collagen substitute did, together with better cytocompatibility, tissue ingrowth, and neoangiogenesis. This composite substitute further demonstrated in clinical case study its ideal repair effect with no CSF leakage or other adverse reactions.

Conclusion: In conclusion, the new biomimetic composite substitute provides alternative substitute for dura repairing.

Keywords: Dural substitute, composite substitute, collagen substitute, tissue regeneration

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
TABLE 1. A great deal of published literature is available which details the biocompatibility of PLA in various applications. Selected articles are organized by various clinical applications in the table below.

Tissue Type	Conclusions	References
PLA and Tissue Reactions	PLA is biocompatible with soft tissues.	3,6,10,12,17,26,43,53,55
PLA and Human Orthopedic Procedures	PLA is biocompatible when used in orthopedic procedures in Humans.	40,41,42,44
PLA and Animal Orthopedic Procedures	PLA is biocompatible when used in orthopedic procedures in Animals.	24,33
PLA and Bone Healing	PLA is biocompatible with bone tissue and does not adversely affect bone healing.	44,49
PLA and Marrow Cells (Osteocytes)	PLA is biocompatible with bone marrow cells.	18,32,34
PLA and Craniofacial Bone Healing	PLA is biocompatible with craniofacial bone tissue and does not adversely affect craniofacial bone healing.	6,29
PLA and Mandibular Bone Healing	PLA is biocompatible with mandibular bone tissue and does not adversely affect mandibular bone healing.	4,49,50
PLA and Periodontal Bone Healing	PLA is biocompatible with periodontal bone tissue and does not adversely affect periodontal bone healing.	8,10,39
PLA and Femoral Bone Healing	PLA is biocompatible with femoral bone tissue and does not adversely affect femoral bone healing.	53,54,55
PLA and Dural Tissue	PLA is biocompatible with dural tissue.	23,30,31,35,57
Spinal Cord and Schwann Cell Biocompatibility	PLA is biocompatible with spinal cord and spinal cord Schwann cells.	16
Effect of Lactic Acid on Brain Tissue, CSF and pH	The cerebrospinal fluid is capable of buffering up to 4 times the anticipated amount of lactic acid with no change in pH. Brain tissue is biocompatible with PLA.	22,25,37
PLA and Nerve Tissues	PLA is biocompatible with peripheral nerves.	37,38,45

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THE USE OF GELFOAM FILM AS A DURAL SUBSTITUTE

PRELIMINARY REPORT

WALTER G. SCHEUERMAN, M.D., FRANCISCO PACHECO, M.D.,
AND ROBERT A. GROFF, M.D.

*Harrison Department of Surgical Research, School of Medicine, University of Pennsylvania,
and Department of Neurosurgery, Graduate Hospital, University of Pennsylvania,
Philadelphia*

From the foregoing experimental and clinical data, it has been shown that gelfoam film is an effective dural substitute. It prevents adhesions from developing between the dura and cortex with practically no tissue reaction. The fact that the film is absorbed in about 80 to 90 days overcomes the disadvantage to all foreign material. The ease of handling and its low cost permits a wide use as compared to fibrin film, yet it has all the advantages of the latter.

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