neon³[™]

universal OCT spinal stabilization

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The perfect solution for
your daily challenge

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Komponenten

Introduction

Indications and contraindications

neon³™ implants

neon³™ instruments

Surgical technique

■ General notes

■ Implantation of lateral mass screws and long shaft screws

■ Implantation of pedicle screws

■ Implantation of transarticular screws

■ Occipito-cervical stabilization

■ Reduction

■ Connectors

■ Cross links

■ Hooks

■ Navigation of the drill guide with Brainlab

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
The small pictures at the bottom of the page show the step-by-step application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-colored.

Art. Nr./Art. No.

Art. Nr./Art. No.

The neon^{3™} spinal system is a modular system for posterior stabilization.

The system offers maximum variability and allows all standard screw placement techniques. Application extends from the occipitocervical junction to the upper thoracic area.

This surgical technique describes the implant, the instruments, and the steps involved in the application of neon^{3™}. It is not sufficient as the sole basis for a successful application of the neon^{3™}. It is recommended to master the surgical technique with an experienced surgeon. Please observe the instructions for use for neon^{3™} implants, the assembly and disassembly instructions and the instrument processing instructions.

Please observe the neon^{3™} surgical technique WS 1402.

Please observe the "Additional information for handles with torque limitation" on page 79.

The development of the neon^{3™} focused on safety during application. The cannulated screws are positioned and anchored securely and precisely using the drill and guide wires. This can also be performed percutaneously as needed using the trocar system. The result is an optimal screw positioning with the greatest possible soft tissue protection. The drill guide is equipped with a special adapter for use of a CAS (Computer Assisted Surgery) system.

For the C1/C2 fixation which occurs particularly frequently, the neon^{3™} atlas claw and special self-drilling screws are available. Occipitocervical stabilizations can be performed as desired with the occipital plate midline 47 mm, 57 mm or with two separate lateral neon^{3™} occipital plates.

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Intended use

The neon³™ occipito-cervical system is used for surgical stabilization, fixation and correction of misalignment of the cervical spine in humans and includes the occipitocervical and cervicothoracic transitions if necessary.

Indications

Occipitocervical, cervical or cervicothoracic, mono- or multisegmental instabilities of rheumatic, traumatic, neoplastic, congenital and iatrogenic genesis, as well as degenerative diseases.

Contraindications

- Substance defects of the anterior spine without additional anterior support
- Patients with a history of material allergy or who tend to react to foreign bodies

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Rod, titanium alloy,
Ø 4.0 mm, length 240 mm,
with hex end

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Rod, CoCr, Ø 4.0 mm,
length 400 mm, with hex
end, long hole

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Rod, pure titanium,
Ø 4.0 mm, length 240 mm,
with hex end, long hole

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
Lateral mass (LMS), long shaft (LSS), pedicle and transarticular (TAS) screws

The different polyaxial screw types (LMS, PS and TAS) are identified by different screw head colors. The color-coding of the lower parts of the screws clearly indicates the diameter of the screw. It is recommended to implant all cannulated screws using a guide wire.

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
Screw, Ø 3.5 mm,
lengths 10 - 20 mm

Properties

- Lateral mass application, C1 to C7
- Polyaxial movement $\pm 48^\circ$
- Not cannulated
- Self-tapping

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
Screw, Ø 4.0 mm,
lengths 10 - 20 mm

Properties

- Lateral mass application, C1 to C7
- Polyaxial movement $\pm 46^\circ$
- Not cannulated
- Self-tapping

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CS 3904-40-xx


Screw, long shaft, cannulated, Ø 4.0 mm,
lengths 26 - 40 mm

Properties

- Lateral mass application, C1
- Polyaxial movement $\pm 46^\circ$
- Cannulated
- Self-tapping

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
Screw, cannulated, Ø 4.0 mm,
lengths 20 - 36 mm

Properties

- Transpedicular, cervical application (C2+C7)
- Polyaxial movement $\pm 38^\circ$
- Cannulated
- Self-tapping

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CS 3906-40-xx

Screw, cannulated, favored angle medial, Ø 4.0 mm,
lengths 20 - 36 mm

Properties

- Transpedicular, cervical (C3-C6), translaminar (C2) application
- Polyaxial movement $-11^\circ/+55^\circ$
- Cannulated
- Self-tapping

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
Screw, cannulated, Ø 4.5 mm,
lengths 25 - 50 mm

Properties

- Transpedicular, thoracic application (Th1-Th3)
- Polyaxial movement $\pm 30^\circ$
- Cannulated
- Self-tapping

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CS 3907-55-xx

Screw, cannulated, Ø 5.5 mm,
lengths 25 - 55 mm, perforated

Properties

- Transpedicular, thoracic application (Th1-Th3)
- Polyaxial movement $\pm 30^\circ$
- Cannulated and perforated
- Self-tapping

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CS 3909-40-xx

Screw, cannulated, favored angle cranial, Ø 4.0 mm, lengths 34 - 50 mm

Properties

- Transarticular C1/C2 application
- Polyaxial movement +55°/-11°
- Cannulated
- Self-drilling

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Occipital plate, midline,
47 mm, for OC rod

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Occipital plate, midline,
57 mm, for OC rod

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OC rod, titanium alloy,
Ø 4.0 mm, straight,
length 240 mm, long hole

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OC rod, titanium alloy,
Ø 4.0 mm, straight,
length 240 mm,
screw hole

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Occipital screw, Ø 5.0 mm,
lengths 6 - 16 mm,
self-tapping

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Counter hook, long,
for atlas claw, Ø 4.0 mm

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Offset connector, parallel,
right, lengths 7 - 20 mm

CS 3918-07-xx

Offsetverbinder, parallel,
links, Längen 7 - 20 mm

Offset connector, parallel,
left, lengths 7 - 20 mm

CS 3919-01

Stab-Stabverbinder,
zervikothorakal, axial

Rod-rod-connector,
cervicothoracic, axial

CS 3919-02

Stab-Stabverbinder,
zervikothorakal, parallel

Rod-rod-connector,
cervicothoracic, parallel

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***Note:** The instrumentation of these implants is performed analogously to the neonTM connectors (CS 1409-15T, -17T, -19T), whose application is described in the surgical technique neonTM WS 1402. CS 3918-05-01, -02, -03 must be stored in the implant tray of the neonTM system, however they must be separate from the neonTM connectors.

Warning:

Caution: Risk of mix-ups! The neonTM connectors CS 1409-15T, -17T and -19T are **NOT** suitable for use with rods from the neon^{3TM} system (basic rod Ø 4.0 mm). The neonTM connectors CS 3918-05-01, -02 and -03 are **NOT** suitable for use with rods from the neonTM system (basic rod Ø 4.5 mm).

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****Note:** The instrumentation of these implants is described in neonTM surgical technique WS 1402.

*****Note:** The locking washer is used as an anti-twist feature in the OC plate of neon^{3TM} (alternatively to CS 1417-2T in neonTM).

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Probe, cannulated, for
screws Ø 4.5 / 5.5 mm,
coupling B

CS 3921-02

Trokardraht
für Ahle, thorakal

Trocar tipped wire
for probe, thoracic

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Trocar tipped drill wire,
Ø 1.5 mm, length 410 mm

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Trocar, cannulated, for
trocar drill wire, Ø 1.5 mm

CS 3923-02

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Sleeve for trocar tipped
drill wire

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Tap, cannulated,
for screws Ø 5.5 mm,
coupling B

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Reference star for
use with Brainlab

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Screwdriver with SR 20,
self-retaining

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Introducer for
rod-rod-connector

CS 3950-01

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OC drill guide for drill depth
6 and 8 mm

CS 3950-02

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OC drill guide for drill depth
10 and 12 mm

CS 3950-03

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OC drill guide for drill depth
14 and 16 mm

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OC drill, flexible,
coupling AO
**Caution:
For single use only!**

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OC screwdriver,
flexible, self-retaining,
coupling AO

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Adapter for injection
(dead volume 0.3 ml)
**Caution:
For single use only!**

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Handle S with torque limitation, 3,5 Nm, coupling AO

Note: Please observe the "Additional information for handles with torque limitation" on page 79.

Note: Please observe arrow markings when used in combination with CS 3926-03, CS 3927-03, CS 3944-02, CS 3951-01/-02.

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Pedicle probe with ball handle, thoracic, length 230 mm

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Handling the instrument set and the introduction of the implants are shown using a model of the spine.

The general and recurrent steps are summarized in the following chapter.

- Assembly of the drill guide and trocar
- Depth determination
- Screw implantation
- Bending the rods
- Final fixation

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Assembly of the drill guide

The actuator of the drill guide (CS 3929-05) is turned clockwise [1]. The depth stop (CS 3929-02) can now be pushed into the drill guide (CS 3929-05) [2].

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Adjusting the drill depth

By pushing the actuator, the depth stop (CS 3929-02) can be adjusted in the drill guide to the desired depth and can be secured by releasing the actuator. Depth adjustment possible in 2-mm intervals.

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

Attaching the trocars and adapters for drill on the drill guide

The trocars (CS 3923-01, CS 3928-04), adapters for drill (CS 3926-04, CS 3926-05, CS 3927-04) and drills (CS 3926-03, CS 3927-03) are color-coded and are oriented to the screw head. Adapters for drills and trocars can be attached using a coupling on the drill guide (CS 3929-05).

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Depth measurement with the depth gauge

The depth gauge (CS 3930) is used to check the drill depth.

In the case of bicortical screw placement, the opposing cortical bone can be felt using the probe hook. The depth can now be read on the instrument.

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Note: When long shaft screws (CS 3904-40-xx) are used, it should be borne in mind that the threadless portion of the lower part of the screw is always 13 mm long and this should be subtracted from the total length indicated in order to determine the correct screw variant. See also figure on page 8.

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
Instruments for screw implantation

When the self-retaining protection sleeve (CS 3932-02) is needed, this sleeve must be pushed onto the screwdriver (CS 3932-04) before securing the screw. Then the screw is secured on the instrument.

The screwdriver (CS 3932-04) is to be used for all screw types of the neon³™ system. It is cannulated and can be pushed over a drill wire or guide wire (CS 3922-410, CS 3924-410).

If the screwdriver is be used similarly to a watchmaker's handle, the rotating end, self-retaining (CS 3932-05) on the end of the handle can be inserted in the handle.

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
Aligning the screw heads

The heads of the inserted screws can be aligned using the screwhead adjuster (CS 3933) that fits the length of the rod.

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Determining rod length

The rod length can be determined using the measuring calipers (CS 5788). Alternatively, the length and shape of the rod can be determined using the phantom rod (CS 3934).

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Inserting the phantom rod

The phantom rod (CS 3934) can be inserted into the aligned screw heads using the rod holder (CS 3940-01) and bent into shape.

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Cutting the rods

The neon³™ rods made of a titanium alloy are available in graduated lengths. If these rods or rods made of CoCr, OC rods, the rod for the atlas claw or the pure titanium rod need to be shortened, cutting pliers must be used (not included in the scope of the system).

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Bending the rods

The neon³™ rods can be shaped appropriately using the rod bending forceps (CS 3936-01). In the basic setting (button A and B are both down ↓), small angles of bending can be made along the rod.

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If bending angles up to 90° are needed, the following procedure is to be used (button down ↓, button up ↑):

	A	B
Step 1	↓	↓
Step 2	↑	↑
Step 3	↑	↓
Step 4	↓	↓

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Note: The individual bending steps must be completely performed to enable rod insertion into the instrument for the following step.

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Placement of the locking screw

The self-retaining screwdriver (CS 3944-03) is used for the pre-fixation of the locking screw (CS 3901-01) on the screw head.

Note: Before screwing in, the locking screw must be rotated a half-turn in the opposite direction to ensure alignment of the locking screw to the screw head.

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Final tightening of the locking screw

S-handle with torque limitation, 3.5 Nm (CS 7154) can only be used with the non-self-retaining screwdriver with SR 20 (CS 3944-02) with AO coupling.

The counterpart (CS 3945-01) is inserted up to the stop on the screw-rod connection. Then the locking screw (CS 3901-01) can undergo final tightening with the screwdriver (CS 3944-02), which is guided through its guide barrels centrally in the counterpart (CS 3945-01). The handle S with torque limitation (CS 7154) releases at a torque of 3.5 Nm.

Note: Please observe the "Additional information for handle with torque limitation" on page 79.

■ Implantation von Massa lateralis- und Langschaftschrauben ■ Implantation of lateral mass screws and long shaft screws

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1. Screw placement technique

Screw placement in the lateral mass (C3 to C7)

Screw placement in the lateral mass C3 to C7 is the standard instrumentation technique due to the simple and secure application.

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Magerl technique


The entry point of the screw is 1-2 mm medially and cranially from the center of the lateral mass. The screw should be inserted 20 - 30° diverging laterally and parallel to the upper joint facet in the direction of the anterolateral border of the upper joint process. In this way, a screw length as large as possible with a low risk of injury to neural and vascular structures is possible.

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
An's technique

The point of entry of the screw is in the middle of the lateral mass. The screw should be inserted 30° diverging laterally and at a 15° upward gradient.

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
Screw placement in the lateral mass in C1, in the isthmus in C2

In the case of screw placement in the lateral mass in C1, the screw runs in the sagittal plane at a slight upward gradient in relation to the C1/C2 joint gap and at a right angle to the frontal plane. The point of entry is beneath the dorsal atlas arch at the junction of the same to the lateral mass C1. The C2 root, which runs posterior to the point of entry, is to be held at a distance in a caudal direction using a dissector prior to drilling and the introduction of the screw.

The isthmus screw has the same point of entry and the same pathway in C2 as the transarticular screw C1/C2. However, it ends at the cranial joint surface of C2 or at the foramen of the vertebral artery, provided there is a "high riding transverse foramen."

- Implantation von Massa lateralis- und Langschaftschrauben
- Implantation of lateral mass screws and long shaft screws

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2. Surgical technique for lateral mass and long shaft screws

Opening the cortical bone

After identifying the point of entry, this point is decoricated using the awl (CS 3920-01) or a high-speed cutter.

Preparing the screw hole

As described under 5.1., the drill guide (CS 3929-05) is mounted with the adapter for drill 1 (CS 3926-04 or -05) and drill 1 (CS 3926-03). For better identification, the drill and adapter for drill are marked with a **green ring**.

The drill guide is then used to prepare the screw hole for screws CS 3902-xx-xx, CS 3903-xx-xx and CS 3904-xx-xx.

Tapping (optional)


As support for pedicle preparation in the case of sclerotic cancellous vertebral bone, the tap (CS 3925-02 for screws Ø 3.5 mm, CS 3925-03 for screws Ø 4.0 mm) with the S-handle (CS 7152) can be used.

The further procedure is then performed analogously to the steps described in the general notes (chapter 5.1.).

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
1. Transpedicular screw in the cervical spine

Instrumentation with pedicle screws leads to significantly increased primary stability of the instrumentation in comparison to the instrumentation with lateral mass screws.

To determine the screw entry points, computed tomography with multiplanar reconstructions of the pedicle of the vertebra to be instrumented should be performed preoperatively. Using multiplanar reconstructions, the point of entry, the angle of inclination and the length of the screw can be preoperatively determined. In general, the point of entry from C3 to C7 is just below the caudal joint facet of the cranial vertebra, far laterally.

The convergence of the screws is 40 to 50° on average but may considerably deviate in individual cases. The screw direction in C2 deviates somewhat with a small amount of convergence. In the case of a short-segment instrumentation and small access, percutaneous screw insertion is useful for achieving the necessary convergence.

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Implantation of pedicle screws without a sleeve (open application)

Opening the cortical bone


After identifying the point of entry, this point is decorated using the awl (CS 3920-01) or a high-speed cutter.

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Mounting of the drill guide

As described under 5.1., the drill guide (CS 3929-05) is mounted with the adapter for drill 2, serrated (CS 3927-04) and mounted to drill 2 (CS 3927-03). For better identification, drills and adapter for drill are marked with a **brown ring**.

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
Implantation of pre-angled pedicle screws with trocar system and sleeve (percutaneous application)

Opening the cortical bone

After identifying the point of entry, this point is decor-ticated using the awl (CS 3920-01) or a high-speed cutter.

Mounting of the drill guide


As described under 5.1., the drill guide (CS 3929-05) is mounted with the trocar for drill 2 (CS 3928-04) and drill 2 (CS 3927-03). In addition, the sleeve for drill 2 (CS 3928-05) still needs to be guided over the trocar. For better identification, drills and trocar are marked with a **brown ring**.

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Preparation of the pedicle

For the transpedicular screws in the area of C3 to C6, the drill hole is made using the trocar system (CS 3928-04, -05).

The trocar system is inserted through lateral stab incisions and drill 2 (CS 3927-03) is introduced. After measuring the screw length using the depth stop, the drill is removed and the blunt guide wires (CS 3924-410) are introduced through the trocar. The length of the screws can be measured once more, as needed, by sliding the depth stop up to the laser marking on the blunt guide wire. Then the drill guide (CS 3929-05) is removed and in the process, the guide wires and sleeves remain in situ.

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Placement of the favored angle pedicle screws


As already described in the general instructions, the screw is secured on the screwdriver (CS 3932-04). The direction of the favored angle is shown by the laser markings on the screw head [photo]. Then the screw is inserted via the guide wire (CS 3924-410) through the sleeve for drill 2 (CS 3928-05).

Thereafter the screw heads are aligned, the rod appropriate for the length is inserted, and the locking screws (CS 3901-01) are inserted using the mounted handle S with torque limitation, 3.5 Nm (CS 7154 and CS 3944-02) and counterpart (CS 3945-01) as described in the general notes and finally secured with 3.5 Nm.

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2. Translaminar screw placement in C2

Wright's technique

The translaminar screw placement of C2 with cross-over screws placed bilaterally is a rigid and technically simple fixation.

The point of entry is at the junction of the spinous process of C2 to the right lamina, near the anterior border of the lamina of C2.

The screw placement in the opposite (left) lamina is parallel to the exposed opposite surface of the lamina whereby the angle is selected such that a possible cortical perforation preferably occurs on the posterior side.

The second screw is placed in the right lamina of C2 with the screw entry point at the junction of the spinous process of C2 to the left lamina, along the caudal border of the lamina.

Note: Due to the special screw head orientation the use of the flexible screw driver CS 3955-01 might be indicated, which is stored on layer 1 for neon^{3™} OC-module.

Note: The use of offset connectors, parallel is described on page 58.

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
3. Stable-angle pedicle screws in the cervical spine

Note: Based on the nominal breaking point on Easyfit, red (CS 1413), for use without a spacer, the use of spacers (CS 1412-2T and CS 1412-4T) is not permitted.

Note: The instrumentation of these implants is performed analogously to the neonTM connectors (CS 1409-15T, -17T, -19T), whose application is described in the surgical technique neonTM WS 1402. CS 3918-05-01, -02, -03 must be stored in the implant tray of the neonTM system.

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4. Transpedicular screws in the thoracic spine

Transpedicular screw placement in the thoracic spine

In the upper and mid thoracic spine, the point of entry of the screws is generally on the upper edge of the transverse process, at the deepest point of the concavity. The convergence is 0 to 15° in most cases. The direction of the screws slightly ascends ventrally in relation to the cover plate of the vertebra.

Preoperative planning of the screws using computed tomography with multiplanar reconstructions is recommended.

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
Assembling the cannulated awl

The gearshift handle (CS 7150) is positioned on the cannulated awl (CS 3921-01) and the trocar wire (CS 3921-02) is screwed in. The same awl is used to prepare both thoracic screws with diameters of 4.5 mm and 5.5 mm.

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
Preparation of the pedicle

The cortical bone of the pedicle is opened at the entry point using the pre-assembled awl. The screw hole is prepared further, whereby the awl is driven further through the pedicle into the vertebral body using constant rotation until it reaches the later intended position of the tip of the screw. Lateral bending of the awl should be avoided. The trocar wire is subsequently unscrewed and removed and the awl remains in the pedicle. The awl can likewise be used to determine screw length.

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Inserting the guide wire

The blunt guide wire (CS 3924-410) is introduced through the cannulation of the awl. The awl is then carefully removed while the guide wire remains in the pedicle.

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
Tapping (optional)

To support pedicle preparation in the case of sclerotic cancellous vertebral bone, the tap (CS 3925-04 for screws Ø 4.5 mm or CS 3925-05 for screws Ø 5.5 mm, cannulated) can be used.

Important: When using the tap CS 3925-05, it must be held in position during application via the guide wire so that it is not pushed forward. Lateral C-arm control.

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Placement of the thoracic pedicle screws

The thoracic pedicle screws (CS 3907-45-xx and CS 3907-55-xx) are inserted using the screwdriver (CS 3932-04). For details on handling the screwdriver (CS 3932-04), see general notes.

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Cement augmentation

In the neon³™ system, the thoracic pedicle screws with Ø 5.5 mm (CS 3907-55-xx) are provided with lateral perforations for cement augmentation.

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Mounting the adapter

Cement augmentation via the perforated pedicle screw requires that a guide wire be used. It is used to facilitate assembly and ensures the patency of the screw cannulation.

The adapter (CS 3962-01) is placed on the screw via the guide wire (CS 3924-410) and secured in the screw head. Then the guide wire can be removed. All adapters should be secured to the screws before starting the augmentation.

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Applying the cement

Mix the cement according to the instructions and draw it into the syringes. The syringe is secured to the adapter via the Luer lock end piece. Ideally, all syringes are positioned on the screws to be augmented.

Apply the cement after the waiting time has elapsed and the cement has reached the correct viscosity. The cement may need to be applied quickly depending on how fast it dries.


Important: Use a C-arm. The adapter is intended only for one-time use. To prevent cement back flow, leave the adapters on the screws until the cement has dried.

Then the adapter is released using the handle for injection adapters (CS 3850-03) and removed.

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1. Screw placement technique

Transarticular screw placement in C1/C2

A computed tomography with multiplanar reconstructions should be performed preoperatively before each transarticular instrumentation. Using the CT, a decision should be made as to whether a secure transarticular screw placement C1/C2 is possible, in particular without injury to the vertebral artery.

The point of entry for transarticular screwing C1/C2 is on the straight, sagittal line which leads centrally through the isthmus (connection between the cranial and caudal joint fact) on the lower edge of the caudal joint process C2, 2 mm cranially and laterally of the medial edge.

■ Implantation von transartikulären Schrauben (TAS) ■ Implantation of transarticular screws (TAS)

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2. Surgical technique for transarticular screws (TAS)

The surgical technique is described using the example of atlantoaxial screw placement of C1/C2 according to Magerl in connection with an atlas claw.

Mounting of the drill guide

The trocar, cannulated, for trocar drill wire Ø 1.5 mm (CS 3923-01), the trocar tipped drill wire Ø 1.5 mm (CS 3922-410) and the sleeve for trocar tipped drill wire (CS 3923-02) are used to prepare the screws CS 3909-xx-xx. The sleeve can be secured after attaching the trocar in a self-retaining fashion on the drill guide (CS 3929-05). The drill wire has a circumferential laser marking which, in connection with the depth stop (CS 3929-02), shows the drill depth. The drill depth corresponds to the screw length.

Beidseitig eingebrachte Bohrdrähte (CS 3922-410)
Parallel liegende Trokarspitzen vor der Schraubeneinbringung
Trocar tipped drill wires (CS 3922-410) introduced bilaterally
Trocar sleeves positioned in parallel prior to screw introduction

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Preparation of the screw hole

To keep the access small in the case of short-segment instrumentations, the percutaneous instrumentation with the trocar system via two upper thoracic additional incisions is recommended. The trocar system, consisting of a sleeve for trocar tipped drill wire (CS 3923-02), trocar, cannulated, for trocar drill wire Ø 1.5 mm (CS 3923-01), is mounted on the drill guide and inserted via a skin incision. The trocar tipped drill wire, Ø 1.5 mm (CS 3922-410) is introduced into the drill guide and applied transarticularly. The length of the screw can be read on the depth stop (CS 3929-02).

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Atlantoaxial instrumentation of a transarticular screw (TAS)

Opening the cortical bone

If necessary, the cortical bone can be opened using the cutter with handle, cannulated, Ø 3.0 mm (CS 3925-01). In doing so, this is pushed through the sleeve for trocar tipped drill wire (CS 3923-02) via the trocar tipped drill wire (CS 3922-410) and the cortical bone is opened.

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Implanting the transarticular screw (CS 3909-40-xx)


The self-drilling transarticular screw has a cranial favored angle. The screw with corresponding length is mounted on the screwdriver, cannulated, with SR 15 (CS 3932-04) and pushed over the trocar tipped drill wire (CS 3922-410) and implanted. The direction of the favored angle is shown by the laser markings on the screw head. Because of the tip of the drill, the screws should in no case be too long in order to avoid irritation to the prevertebral soft tissues. The drill wire is removed after screw implantation.

Note: When implanting the screw, the drill wire must be kept in position. Lateral C-arm control.

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Mounting and inserting the neon³™ atlas claw

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The rod for atlas claw, Ø 4.0 mm (CS 3915-01) is shortened to the required length and the appropriate counter hook for atlas claw (CS 3915-02 or CS 3915-03) is positioned on the end of the rod. The rod for atlas clamp is suspended on the vertebral arch of the atlas, inserted in the screw head of the TAS and prefixed with the locking screw (CS 3901-01) on the screwdriver with SR 20 (CS 3944-03).

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Note: The rod for the atlas claw should not be fixed further than 16 mm laterally from the midline in order to avoid injuries to the vertebral artery.

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Positioning the counter hook for neon³™ atlas claw

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After final tightening the locking screw (CS 3901-01) on the TAS (see General notes page 29), the counter hook (CS 3915-02 or -03) is pushed with the compression forceps (CS 3943-04) under the caudal end of the atlas vertebral arch and the atlas claw is closed.

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The compression forceps (CS 3943-04) are supported on the screw head of the TAS and the counter hook (CS 3915-02). The locking screw of the counter hook is tightened hand-tight using the screwdriver with SR 20 (CS 3944-02) and the handle S with torque limitation (CS 7154).

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The occipitocervical (OC) stabilization can be performed with three different implants, the occipital plate 47 mm (CS 3913-06), 57 mm (CS 3913-07), or two lateral (CS 3913-05) neon^{3™} occipital plates. The neon^{3™} occipital plates, lateral, are easier to mount, particularly in the case of anatomically difficult situations.

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1. Instrumentation of medial occipital plates

Check plate position and the required bending of the plate

The plate is held on the occiput in order to determine the required bending and position of the plate.

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Inserting the plates in the plate bender

The plate can now be inserted in the plate bending clamp (CS 3954-01). To secure the plate to the plate bending clamp, the actuator flap of the plate bending clamp must be pushed forward. With the actuator pressed, the plate can then be fully pushed into the instrument. By releasing the flap, the plate is secured to the instrument.

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Bending the lateral flaps

The lateral flap of the plate is inserted all the way into the larger recess on the OC plate bender with handle (CS 3954-03) (see photo, left). The plate can now be bent into the desired position.

Warning: Repeated bending back and forth must be avoided.

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
Bending the upper screw flap

The upper screw flap can be bent like the lateral flaps. To avoid pressure marks or damage to the back-out mechanism, the smaller recess on the plate bender (CS 3954-03) must be used for bending the upper screw flap (see photo, left).

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
Pre-fixation of the plate

The pre-bent plate is positioned in situ and pre-fixed on the occiput using 2 OC pins for plate fixation (CS 3957-01). For this purpose, the pin is connected with the screwdriver with SR 20, self-retaining (CS 3944-03) and secured to the occiput through the screw hole of the plate.

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Preparation of the screw hole on the occiput

To drill the screw holes on the occiput, a rigid (CS 3951-01) or a flexible drill (CS 3951-02) can be selected. When drilling by hand, the respective drill must be connected with the handle S with torque limitation (CS 7154). The OC drill guide (CS 3950-01, -02, -03) is divided into 2-mm increments from 6-16 mm. Through a spherical surface on the OC drill guide, the screw holes can be created a maximum of $\pm 8^\circ$ obliquely.

Fixed drill depth of:
6 mm + 8 mm (CS 3950-01)
10 mm + 12 mm (CS 3950-02)
14 mm + 16 mm (CS 3950-03)

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
Screwing the medial plate:

The occipital screws (CS 3914-50-xx) can now be inserted with the aid of the screwdriver with SR 20, self-retaining (CS 3944-03) or the OC screwdriver, flexible, self-retaining (CS 3955-01 + CS 7154). Then the OC pins for plate fixation are removed once more. The freed-up screw holes can now be drilled, as previously shown, and the remaining occiput screws can be inserted.

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Inserting the TAS screws

The TAS is inserted as described in the general notes.

Attaching the rods to the occipital plate

There are several straight OC rod variants available:

Rods in titanium alloy (TiAl6V4), length 240 mm

CS 3913-40-03 OC rod, titanium alloy, Ø 4.0 mm, straight, long hole
CS 3913-40-04 OC rod, titanium alloy, Ø 4.0 mm, straight, screw hole

Rod in pure titanium

CS 3912-240 Rod, pure titanium, Ø 4.0 mm, length 240 mm with hex end, long hole

Rod in CoCr (optional)

CS 3911-400 Rod, CoCr, Ø 4.0 mm, length 400 mm with hex end, long hole

These must be correspondingly preformed using the rod bending forceps (CS 3936-01).

Warning: Bending reduces the mechanical stability of the instrumentation.

The OC locking screw (CS 3913-04) is placed in situ using the self-retaining screwdriver (CS 3944-03) and pre-fixed. Then, using the screwdriver (CS 3944-02) and the associated torque handle (CS 7154), the screw (CS 3913-04) undergoes final tightening with 3.5 Nm. In the event of difficult space conditions, the flexible screwdriver (CS 3955-01) can be used.

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Instrumentation with lateral OC plates with rod



Lateral occipital plate

The surgical technique for the occipital plate, lateral, with rod (CS 3913-05) differs from that of the plates to be fixed medially (CS 3913-06/-07) only with regard to the bending process.

In the process, the plate is held firmly at the lower caudal screw hole with the OC plate bending clamp (CS 3954-01), as shown in Fig. [1]. The bending process at both bending sites is then performed using the OC plate bender with handle (CS 3954-03) in that the plate is introduced into the smaller recess (see Fig. [2]). The bending sites are suitable for bending in the frontal as well as in the sagittal plane. The rod of the lateral plate must be brought into the desired shape using the rod bending forceps (CS 3936-01).

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1. Reduction with the persuader

Caution: Not recommended for the use of LMS screws

Mounting the rod pusher (persuader)

The inner sleeve for rod pusher (CS 3938-04) is inserted into the outer sleeve for rod pusher (CS 3938-05) such that the pins can be guided into the recess [1]. If the first stop has been reached, a right rotation of the inner sleeve is made and after releasing it, the sleeve snaps securely into place in the outer sleeve.

The mounted rod pusher (CS 3938-04/-05) is suspended on the screw head with the outer sleeve pulled back from below [2]. Then the outer sleeve is released, and this sleeve is thus placed on the rod [3]. By rotating the knurled nut of the inner sleeve for rod pusher (CS 3938-04), the rod is guided completely through the outer sleeve into the screw head [4]. The end position is shown through the markings on the inner and outer sleeve [5].

Caution: A reduction of up to 12 mm can be performed using this instrument.

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Screwing in the locking screw

The locking screw (CS 3901-01) is inserted through the rod pusher (CS 3938-04, -05) using the screwdriver with SR 20, self-retaining (CS 3944-03) and pre-fixed.

The rod pusher can be removed by re-releasing the outer sleeve from the screw head.

Final tightening with the counterpart (CS 3945-01) and the handle S with torque limitation, 3.5 Nm (CS 7154), as described in the general notes.

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
2. Reduction via the rocker

Caution: Not recommended for the use of LMS screws

Insertion of the rod with the rod pusher (rocker)

The fork of the rod pusher (rocker) (CS 3939-01) is suspended on the screw head from below. By pivoting the rocker, the rod can be inserted into the screw head.

Caution: A reduction of up to 5 mm can be performed using this instrument.

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Placement of the locking screw

The locking screw (CS 3901-01) is inserted using the screwdriver with SR 20, self-retaining (CS 3944-03).

Final tightening with the counterpart (CS 3945-01) and the handle S with torque limitation, 3.5 Nm (CS 7154), as described in the general notes.

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1. Offset connector, straight and parallel

Two straight and twelve parallel variants of offset connectors are available in the event of an overly wide offset between the rod and the screw head.

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Offset connector, square (CS 3918-02-02 or CS 3918-02-03, long)

The screwdriver with SR 20, coupling AO (CS 3944-02) is initially mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). Thereafter, the appropriate offset connector, square, is placed on the rod using the self-retaining screwdriver (CS 3944-03) as well as inserted in the screw head and then pre-fixed. The locking screw CS 3901-01 is now prefixed on the screw head using the self-retaining screwdriver. The grub screw in the offset connector and the locking screw in the screw head are then finally secured with 3.5 Nm using the mounted screwdriver (CS 3944-02 and CS 7154). The rod slot of the offset connector is always aligned at a right angle to the screw head.

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Größenabstufungen: 7 mm
9 mm
11 mm
14 mm
17 mm
20 mm

Offset connector, parallel right, (CS 3918-06-07 to /-20) left, (CS 3918-07-07 to /-20)

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
The screwdriver with SR 20, coupling AO (CS 3944-02) is initially mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). Thereafter, the appropriate offset connector, parallel, is placed on the rod using the self-retaining screwdriver CS 3944-03 as well as inserted in the screw head and then pre-fixed. The locking screw CS 3901-01 is now prefixed on the screw head using the self-retaining screwdriver. The grub screw in the offset connector and the locking screw in the screw head are then finally secured with 3.5 Nm using the mounted screwdriver (CS 3944-02 and CS 7154).

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2. Cable connector of the neon³™ system

Cable connector (CS 3918-01)

The screwdriver with SR 20, coupling AO (CS 3944-02) is initially mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). Thereafter, the cable connector is placed on the rod using the self-retaining screwdriver CS 3944-03 and pre-fixed. The grub screw in the cable connector is then finally secured with 3.5 Nm using the mounted screwdriver (CS 3944-02 and CS 7154).

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3. Parallel rod-rod connectors of the neon^{3™} system

Rod-rod connector, cervicothoracic, parallel

Note: Parallel rod-rod connectors (CS 3919-02) may **only be used in pairs** per rod (see picture). There should be a screw between the connectors. The connector is suitable in both recesses to accept rods with Ø 4.0 mm from the neon^{3™} system and also for rods with Ø 4.5 mm from the precursor neon[™] system.

The rod-rod connector can be held in position with the introducer for rod-rod connectors (CS 3946). The open side of the connector is pushed laterally on the rod that is already in place. The screwdriver with SR 20, coupling AO (CS 3944-02) is initially mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). The locking screws of the rod-rod connector, parallel (CS 3919-02) are secured with the mounted screwdriver with 3.5 Nm.

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4. Axial rod-rod connectors of the neon³™ system

Rod-rod connector, cervicothoracic, axial

The rod-rod connector, cervicothoracic, axial (CS 3919-01) is suitable on both sides to accept rods with Ø 4.0 mm from the neon³™ system and also for rods with Ø 4.5 mm from the precursor neon™ system.


The rod-rod connector can be held in position with the introducer for rod-rod connectors (CS 3946). The screwdriver with SR 20, coupling AO (CS 3944-02) is initially mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). The locking screws (CS 3901-01) are then secured with the mounted screwdriver with 3.5 Nm.

Caution: Each rod end must be connected by 2 locking screws.

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Determining the implant size

The offset of both rods can be determined using the measuring calipers (CS 5788).

Note: The offset of both rods is determined on the external side of the rods.

The size images of the crosslinks (CS 3917-05/-06/-07) on the lid of the box facilitate the selection of the correct implant [1].

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Implantation of the crosslink

The crosslink selected (CS 3917-05,-06 or -07) is positioned from the outside on the two rods of the instrumentation and pressed together.

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
Fixation of the crosslink

The screwdriver with SR 20, coupling AO (CS 3944-02) is initially mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). The locking screws of the crosslink are then secured with the mounted screwdriver with 3.5 Nm.

Note: It should be ensured that the crosslink is fully seated on the rod. The screwdriver should always be kept at a 90° angle to the rod. To prevent tilting, first the outer and then the middle locking screw should be secured.

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Use of hooks

All hook types (CS 3916-01, -02, -03, -04) can be inserted with the hook inserting forceps (CS 3961). Using the inserter for hooks (CS 1466), the positioning of the hook can be prepared in situ ahead of time.

The locking screw (CS 3901-01) is pre-fixed on the hook using the screwdriver with SR 20, self-retaining (CS 3944-03) with the rod already in place. Then the screwdriver with SR 20, coupling AO (CS 3944-02) is mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). Using the mounted screwdriver, the locking screw with 3.5 Nm then undergoes final securing.

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Attaching the reference star for use with Brainlab

The button on the reference star for use with Brainlab (CS 3929-06) must be held down in order to be able to mount the reference star on the drill guide (CS 3929-05). The reference star must be fully mounted on the adapter on the drill guide. The reference star must be mounted in a properly oriented fashion so that the pin on the drill guide engages in the corresponding slots of the reference star and the locking mechanism of the button snaps into place in the corresponding groove of the drill guide. Then the depth stop (CS 3929-02) can be pushed into the drill guide.

After the Brainlab Disposable Reflective Marker Spheres have been screwed on, the instrument is ready for manual calibration. Using the Brainlab Instrument Calibration Matrix (ICM4) and the associated insert for the Brainlab ICM4 (CS 3929-04), the instrument can be easily calibrated. The reference star must be mounted and oriented relative to the drill guide handle such that during use, the navigation camera has a direct and unimpacted view of the reference star.

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The trocars (CS 3923-01, CS 3928-04), adapters for drill (CS 3926-04, CS 3926-05, CS 3927-04) and drills (CS 3926-03, CS 3927-03) are color-coded and are oriented to the screw head. Adapters for drills and trocars can be attached using a coupling on the drill guide (CS 3929-01).

Bohreinsatz/Trokar	Bohrer	Art.-Nr.
CS 3923-01	Bohrdraht, Ø 1,5 mm	CS 3922-410
CS 3926-04	Bohrer 1, Ø 2,3 mm	CS 3926-03
CS 3926-05	Bohrer 1, Ø 2,3 mm	CS 3926-03
CS 3927-04	Bohrer 2, Ø 2,6 mm	CS 3927-03
CS 3928-04	Bohrer 2, Ø 2,6 mm	CS 3927-03

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Diameter of the drill that fits the adapter for drill/trocar in each case

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Preparing the calibration

The trocars and adapters for drills must be calibrated using the Brainlab Instrument Calibration Matrix (ICM4) and Brainlab navigation software.

For calibration of the trocars the insert for the Brainlab ICM4 (CS 3929-04) is required and must be inserted fully without jamming until it snaps into place in the 30-mm hole in the ICM4.

The adapters for drills are calibrated in the adjacent drill hole of the ICM4 for Ø 5 mm.

For the calibration, the trocars and adapters for drills must be inserted fully up to the stop in the corresponding drill hole of the insert or ICM4.

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Note: Detailed information on handling the Brainlab Disposable Reflective Marker Spheres, the ICM4, manual calibration and navigation of the instruments in combination with the navigation software can be found in the corresponding Brainlab user manual.

Warning: The navigated instruments are very precise and very sensitive medical instruments and must be used with extreme care. If the instrument should fall down or be damaged, the calibration must be carefully checked or the instrument must be sent back for inspection. Improper handling can cause serious injuries to the patient.

Warning: Plan the OR system setup and the position of the reference star on instruments and patients prior to the surgery. The navigation camera must have an direct and unimpacted view of all reference stars. The navigation system should be set up as described by Brainlab.

Warning: The navigated drill guide may only be used in connection with the associated drills and drill wires. ulrich medical® and Brainlab are not liable for the use of non-specified instruments.

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Warning: Ensure that the drill or drill wire is used with the adapter for the correct drill and with the appropriate color coding.

Warning: After changing of trocar or adapter for drill, the drill guide must once again be calibrated with the Brainlab ICM4.

Warning: The calibration of the trocars may only be performed using the insert for the Brainlab ICM4 (CS 3929-04), the calibration of the adapters for drill may only be performed in the Ø 5 mm drill hole of the Brainlab ICM4.

Warning: Ensure that bending forces on the calibrated instruments are avoided since this deviation can affect accuracy.

Warning: Do not change the orientation of the drill guide during drilling since this may cause the drill to break within the bone.

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Removing the implant


The neon³™ spinal system is removed in the sequence opposite to the sequence of implantation steps.

First the screwdriver with SR 20, coupling AO (CS 3944-02) must be mounted on the handle S with torque limitation, 3.5 Nm (CS 7154). The mounted screwdriver and the counterpart for torque wrench (CS 3945-01) are then used to remove the locking screws.

After removing the rod, the screw is removed using the screwdriver, cannulated, with SR 15 (CS 3932-04).

For emergency cases, the uRevision™ Set CS 7300 explantation system is available.

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

Additional information for handles with torque limitation

Handles with torque limitation (CS 7154) are used in the neon³™. The following additional information should be borne in mind during use:

General

The use of these handles with torque limitation for an application that does not conform to the intended use can lead to damage and may consequently represent a risk to the patient or user. Ensure that the products are only used by appropriately qualified or trained staff. These products are handled and used without direct control by the manufacturer and remain the responsibility of the user. The user bears exclusive responsibility if any damage should occur due to improper use. These instruments are not measurement devices.

Intended use for handles with torque limitation


The handles with torque limitation are used to limit torque to a factory-adjusted fixed value. Deviation of the actual torque from the nominal torque is $\pm 10\%$ when the instrument is delivered. The torque indicated is valid only for clockwise rotation. When rotating counterclockwise, the torque is much higher and not specified. Use with active or motorized class IIa drive systems (for example, drilling machines) is not permitted. Only manual operation is allowed. The permissible operating temperature of the handles with torque limitation is between 10°C and 40°C.

Application

The handles with torque limitation cannot be disassembled and are delivered in a non-sterile state. It is recommended that the handles with torque limitation be checked with regard to functioning after cleaning and sterilizing and prior to the surgical intervention. In particular, the torque should be checked. The handles with torque limitation may only be used with the tool provided for this purpose. As soon as the torque limit has been reached upon clockwise rotation, the mechanism disengages with an audible “click.” The maximum allowable axial load is 100 N.

Calibration

The nominal torque is guaranteed for three years from the time of delivery, if the instrument is used as intended. To ensure safe and proper use, the operator must promptly return the handle with torque limitation to the manufacturer at the end of this time period for calibration.

WS 3901 R1/2014-05