uCentum™

comprehensive posterior system

Offen

Open

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+ Universal pedicle screw system
for all applications

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


Spinal systems by ulrich medical® stand for quality “Made in Germany.”

They are the successful result of systematic development activities and many years of experience in medical technology. For more than 100 years, we have done our best every day for our customers and for patient health with our proprietary and innovative products.

Leading-edge technology and competence from one source

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Komponenten

Introduction

Indications and contraindications

uCentum™ implants

uCentum™ instruments

Surgical technique

■ Preparation

■ Instrumentation with polyaxial and monoaxial screws

■ Instrumentation with Pre-Fixation screws

■ Instrumentation with Pre-Fixation Long-Arm screws

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The small pictures at the bottom of the page show the step-by-step application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-colored.

Art. Nr./Art. No.

Art. Nr./Art. No.

This surgical technique describes the implant, the instruments, and the steps involved in the open-approach application of the uCentum™ system. It is not sufficient as the sole basis for a successful application of the system. It is recommended to master the surgical technique with an experienced surgeon. Please observe the instructions for use, the assembly and disassembly instructions and the instrument processing instructions at www.ulrichmedical.com.

Each uCentum™ screw is cannulated and thus permits secure and precise placement via guide wires as well as minimally invasive and percutaneous application. The perforations at the tip of the screw enable cement augmentation through the already implanted screw (except for screws with Ø 4.5 mm).

The uCentum™ system is distinguished by a slim, low-profile screw head. It thus allows space for bone grafting and a good view for fusion monitoring.

The same instruments are used for polyaxial, monoaxial and pre-fixation screws. All screw variations can be freely combined.

The polyaxial design of the pre-fixation screws can be blocked by tightening the external locking screw with torque without, however, inhibiting any desired movement that is axial or radial to the rod. In this way, there are numerous options for reduction the spine, e.g. alignment following fractures or derotation of scoliosis.

The monoaxial uCentum™ screw has a higher degree of angular stability as compared to the polyaxial screws. This increased stability allows shorter instrumentation procedures and simplified reduction.

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Intended use

uCentum™ is an implant system with polyaxial, monoaxial, pre-fixation and long-arm screws, straight and preformed rods and cross stabilizers for surgical stabilization, anchorage and correction of the thoracic and lumbar spine.

Indications

Degenerative disk disease, deformity, spinal tumor, spondylolisthesis, pseudarthrosis, spinal stenosis, destruction and instability

Contraindications

- Patients with fever or leukocytosis in the case of non-spine associated infections
- Patients with a documented history of material allergy or who tend to react to foreign bodies
- The physician must consider carefully before treating patients who are in a generally unfavorable medical or psychological state which could be made worse by the procedure.
- Pregnancy

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
Locking screw for screw, pre-fixation

Stäbe | Rods

CS 3807-xxx

Stab, Ø 6,0 mm, gerade
Längen 30 mm – 180 mm,
Standardzugang

Rod, Ø 6.0 mm, straight
lengths 30 mm – 180 mm,
standard approach

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The various screw types as well as their diameters in the uCentum™ implant system are differentiated by color-coding of the screws. The screw types are differentiated by the positioning of the color-coding; the diameters of the screws are differentiated via the color.

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Caution: It should be ensured that the correct locking screw (CS 3801-01 or CS 3801-02) is used.

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All pedicle screws are cannulated and it is recommended to implant them using a guide wire.

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CS 3801-01

Application options:

- Open and percutaneous application
- Variable-angle reduction along the rod

3

Ø 4,5 mm
Ø 4.5 mm

Ø 5,5 mm
Ø 5.5 mm

Ø 6,5 mm
Ø 6.5 mm

Ø 7,5 mm
Ø 7.5 mm

Ø 8,5 mm
Ø 8.5 mm

Ø 10 mm
Ø 10 mm


Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Schraube, monoaxial | Screw, monoaxial

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






CS 3804-xxx-xxx

Position of the color-coding:
Clamping plate

Associated locking screw:
CS 3801-02

Application options:

- Open and percutaneous application
- Stable-angle reduction along the rod

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
CS 3805-xxx-xxx

Position of the color-coding:
Clamping plate

Associated locking screw:
CS 3801-02

Application options:

- Open application
- Stable-angle reduction along the rod

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Awl, cannulated, for screws
In four diameters

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Awl, cannulated,
for screws Ø 10.0 mm,
coupling A

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
Pedicle probe with ball handle, lumbar, straight, length 230 mm

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Pedicle probe with ball handle, thoracic, length 230 mm

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Pedicle probe with ball handle, lumbar, curved, length 230 mm

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Tap, cannulated, for screws Ø 5.5 – 10.0 mm coupling A

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Dilator S and gauge for screw length

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Threaded element for
guidewire holder

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Screwdriver, open,
cannulated, hex 3.5 mm,
coupling A

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Protection sleeve for
longarm screw

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
Screwdriver for locking
screw, pre-fixation, self-
retaining

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Screwdriver for locking
screw, self-retaining

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Outer sleeve for
rod pusher, open
(Persuader)

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Inner sleeve for
rod pusher, open
(Persuader)

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Compression forceps,
curved, for rod Ø 6.0 mm

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Distraction forceps,
curved, for rod Ø 6.0 mm

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Handle for distraction and compression forceps

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Blade for distraction and compression forceps, fig. W, X

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Blade for distraction and compression forceps, fig. Y, Z

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Blade for distraction and compression forceps, fig. M, width 47 mm

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Blade for distraction and compression forceps, fig. L, width 90 mm

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Blade for distraction and compression forceps, M-L

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T-handle with torque
limitation, 9 Nm, coupling V

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Screwdriver for locking
screw, coupling V

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Screwdriver for locking
screw, pre-fixation,
coupling V

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T-handle with torque
limitation, 6 Nm,
coupling H

CS 3848

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Screwdriver with SR 20,
coupling H

CS 3850-01

Injektionsadapter, offen
(Totvolumen 0,3 ml)

Adapter for injection, open
(death volume 0.3 ml)

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
Handling the instrument set and the introduction of the implants are shown using a model of the spine.

Exposure of the spine

The appropriate open approach is selected and the posterior spinal region is exposed.


1. Preparation using a cannulated awl and trocar wire

Assembling the cannulated awls

The gearshift handle (CS 7138) is positioned on the cannulated awl for screws (CS 3018-xx, CS 3823-100) and the trocar wire (CS 3821-01) is screwed in. For the various screw diameters 4.5/5.5 mm, 6.5 mm, 7.5 mm, 8.5 mm and 10 mm, there are appropriate cannulated awls available (CS 3018-04, -05, -06, -07, CS 3823-100).

Preparation of the pedicle

The cortical bone of the pedicle is opened at the entry point using the pre-assembled cannulated awl for screws. The screw hole is prepared further whereby the awl is driven further through the pedicle into the vertebral body using constant rotation until it reaches the later intended position of the tip of the screw. Lateral bending of the awl should be avoided. The trocar wire is subsequently removed [2] and the awl remains in the pedicle. The awl can likewise be used to determine screw length [1].

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Inserting the guide wire

The guide wire (CS 2624) is introduced through the cannulation of the awl [1]. The awl is then carefully removed while the guide wire remains in the pedicle [2].

Determining screw length

The gauge for screw length (CS 3820-01) is slid over the guide wire up to the bone. Using the laser marking on the guide wire, the length of the screw can be read on the gauge. The function of the orange marking on the measurement instrument is only needed for percutaneous application.

Taps (optional)

As support for the pedicle preparation in the case of sclerotic cancellous vertebral bone, the tap (CS 3824-xxx) can be used. The length of the tapping instrument area allows preparation of the screw hole up to the opposing cortical bone.

Important: The guide wire must be held in position so that it is not pushed forward. Lateral C-arm control. The tap for screws Ø 4.5 mm cannot be applied with a guide wire.

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2. Preparation with an Awl and blunt Awl

Opening the cortical bone

The cortical bone is opened at the pedicle entry point using the awl (CS 3033 + CS 7138).

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Assembling the cannulated awls

The gearshift handle (CS 7138) is positioned on the cannulated awl for screws (CS 3018-xx, CS 3823-100) and the trocar wire, blunt (CS 3821-02) is screwed in. For the various screw diameters 4.5/5.5 mm, 6.5 mm, 7.5 mm, 8.5 mm and 10 mm, there are appropriate cannulated awls available (CS 3018-04, -05, -06, -07, CS 3823-100).

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
Preparation of the pedicle

The screw hole is prepared further whereby the awl is driven further through the pedicle into the vertebral body using constant rotation until it reaches the later intended position of the tip of the screw. Lateral bending of the awl should be avoided. The trocar wire is subsequently removed [2] and the awl remains in the pedicle. The awl can likewise be used to determine screw length [1].

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Determining screw length

The gauge for screw length (CS 3820-01) is slid over the guide wire up to the bone. Using the laser marking on the guide wire, the length of the screw can be read on the gauge. The function of the orange marking on the gauge is only needed for percutaneous application.

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

3. Preparation with an Awl and Pedicle Probes

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Preparation of the pedicle

For further preparation of the pedicle, additional instruments are available with the pedicle probes (UL 8525-23, UL 8526-23, UL 8527-23).

Palpating the drilling channel

Use then the pedicle probe (UL 8522-27) to check the integrity of the pedicle.

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
Assembly of the screw and screwdriver

The desired handle, gear-shift handle (CS 7138 or CS 7139) or T-handle with ratchet (CS 7140), is mounted on the screwdriver, open (CS 3825-01). The protection sleeve (CS 3825-02) is put over the screwdriver, open (CS 3825-01) and pushed until the sleeve clicks into place.

The screwdriver is secured in the screw head of the screw to be implanted by turning the shaft wheel. The polyaxial movement of the screw is blocked.

By pressing the button on the retractor, see detail [1], the retractor can be moved, pushed onto the shaft wheel of the screwdriver and thus the screwdriver can be secured against unintentional release.

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Implanting the screw

The uCentum™ screw is now implanted using the screwdriver.

Important: The guide wire must be held in position so that it is not pushed forward when the screw is screwed in. Use of a C-arm is recommended (lateral view).

To release the tight connection between the screwdriver and screw head, the detent of the shaft wheel must be lifted by the retractor. For this purpose, the button is pressed and the retractor is pushed posteriorly. Then the connection can be released by turning the shaft wheel.

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The guide wire holder prevents the guide wire from shifting during screw insertion. When correctly used, it is an important safety feature of the system. The use of the guide wire holder is recommended particularly in the case of bicortical screw placement.

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Assembly of the guide wire holder

Important: The guide wire holder can only be used with the gearshift handle (CS 7138 or CS 7139). Use with the T-handle with ratchet (CS 7140) is not possible!

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The threaded element (CS 3822-01) is screwed into one of the two handles as far as it will go. Then the handle for the guide wire holder (CS 3822-02) is screwed into the threaded element; in doing so, the handle should not be screwed all the way into the threaded element.

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The uCentum™ screw is implanted via the guide wire in place using the screwdriver and mounted guide wire holder.


As soon as the tip of the screw comes into contact with bone, the guide wire is connected to the threaded element (CS 3822-01) by screwing in the handle (CS 3822-02) further [1].

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


Important: While the screw is screwed into the pedicle, the handle of the guide wire holder must be held since it should not simultaneously be rotated [2]. In this way, the guide wire is held in position while the screw is screwed in [3].

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For alignment, screw head adjuster (CS 3827) can be used.

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Alternatively, the length and shape of the rod can be determined using the phantom rod (CS 8032-xx).

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Bending the rod

The uCentum™ rod can be brought into the suitable shape using the rod bending forceps according to Lepine (UT 1639-29).

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
Inserting the rod

The uCentum™ rod is inserted using the rod holder (CS 3024). Alternatively, the rod clamp (CS 8038) can be used for this purpose.

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
Connecting screws and rod

The rod pusher, open (Rocker) (CS 3837) provide support during complete introduction of the rod into the head of the screw.

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Please refer to chapter 5.6 if reduction should become necessary.

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Mounting the Sleeve

Prior to insertion of the locking screw the sleeve for locking screw (CS 3836) is positioned on the screw head. This instrument can also support the complete introduction of the rod into the screw head.

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
Screwing in the locking screw

The uCentum™ locking screw (CS 3801-01) is screwed into the screw using the self-retaining screwdriver (CS 3851-01) through the sleeve for locking screw (CS 3836) and it is gently secured. The rod and also the polyaxial head are fixed.

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

For monosegmental distraction and compression, two forceps are available in the standard instrument set (CS 3852-01 compression forceps and CS 3852-02 distraction forceps).

B) Alternative distraction and compression instruments (option for variable distances from 6 to 160 mm)

Selection of blades for distraction and compression

Five different blades are available. The blades can be freely combined. A distraction or compression width from 6 to 160 mm can be achieved.

CS 3032-05 blade, W, X, width 3 or 10 mm


CS 3032-07 blade, Y, Z, width 25 mm

CS 3032-08 blade, M, width 47 mm

CS 3032-10 blade, M-L, width 65 mm

CS 3032-09 blade, L, width 90 mm

For blade CS 3032-05, two different widths can be selected since these can be mounted in both directions [1], [2].

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
For distraction or compression, the handles (CS 3032-00) are inserted in the distraction forceps (CS 3032-01) or compression forceps (CS 3032-02). While inserting or removing the handles, press the button on the instrument (orange arrow in figure).

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
Distraction

For distraction, the single-part (CS 3852-02) or assembled distraction forceps (CS 3032-xx) are used. Distraction is performed gradually in the case of the single-part forceps. By pulling back the toothed rod, the forceps can be once again loosened and removed.

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Compression

For compression, the single-part (CS 3852-01) or assembled compression forceps (CS 3032-xx) are used. Compression is performed gradually in the case of the single-part forceps. By pulling back the toothed rod, the forceps can be once again loosened and removed.

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

Assembly of the torque limiter

The screwdriver for locking screws (CS 3851-02) is inserted into the T-handle with torque limitation (CS 7146). The sleeve of the T handle must be pulled back during insertion (orange arrow in the figure).

Caution: Observe the instructions for use for the specifications for testing the torque of the instruments CS 7146 and CS 7148.

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
Securing the instrumentation

The uCentum™ locking screw is secured with the torque limiter (CS 7146 + CS 3851-02) with 9 Nm in the case of an assembled distraction or compression forceps. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3846 + CS 7144) is used. If the screwdriver is fully seated in the inner grub screw of the locking screw, the colored black ring then shows the correct position of the locking screw in the head of the screw.

Information: Use of a crosslink can be found in chapter 5.5

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Use: In the case of TLIF and PLIF as well as burst fractures and scoliosis, since parallel distraction and compression and rod rotation are possible.

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The protection sleeve (CS 3825-02) is put over the screwdriver, open (CS 3825-01) and pushed until the sleeve clicks into place.

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The desired handle, gear-shift handle (CS 7138 or CS 7139) or T-handle with ratchet (CS 7140), is mounted with the screwdriver, open (CS 3825-01).

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The guide wire holder prevents the guidewire from shifting during screw insertion. When correctly used, it is an important safety feature of the system. The use of the guidewire holder is recommended particularly in the case of bicortical screw placement.

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The threaded element (CS 3822-01) is screwed into the gearshift handle as far as it will go. Then the handle for the guidewire holder (CS 3822-02) is screwed into the threaded element; in doing so, the handle should not be screwed all the way into the threaded element.

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
Screwing in the screw with the guide wire holder (optional)

The uCentum™ screw is implanted via the guidewire in place using the screwdriver and mounted guide wire holder.

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
Alignment of the screw heads

The screw heads of the screws are aligned prior to inserting the rod. For alignment, screw head adjuster CS 3827 can be used.

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Determining rod length

The rod length can be determined using the measuring caliper (CS 5788). Alternatively, the length and shape of the rod can be determined using the phantom rod (CS 8032-XX).

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


Connecting screws and rod

The rod pusher, open (Rocker) (CS 3837) provides support during complete introduction of the rod into the head of the screw.

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Please observe chapter 5.6 if reduction should become necessary.

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

Screwing in the locking screw

Important: For all pre-fixation screws, locking screw for screws, pre-fixation (CS 3801-02), must be used.

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
Assembly of the torque limiter

The screwdriver for locking screw, pre-fixation (CS 3847) is inserted into the T-handle with torque limitation (CS 7146). The sleeve of the T handle must be pulled back during insertion (orange arrow in the figure).

Caution: Observe the instructions for use for the specifications for testing the torque of the instrument CS 7146.

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Assembly of the counterpart


The handle L, M8 (CS 7144) must be screwed onto the counterpart, open (CS 3846).

Clamping the polyaxial motion

The uCentum™ locking screw (CS 3801-02) is secured with the torque limiter (CS 3847 + CS 7146) with 9 Nm. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3846 + CS 7144) is used.

Important: Only the polyaxiality of the head is fixed now. The screw head can still be freely positioned along the rod.

If the screwdriver is fully seated in the inner grub screw of the locking screw, the colored black ring then shows the correct position of the locking screw in the head of the screw.

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Parallel distraction and compression

A) Standard instruments for distraction and compression

For monosegmental distraction and compression, two forceps are available in the standard instrument set (compression forceps CS 3852-01 and distraction forceps CS 3852-02).

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Assembly of the torque limiter

The screwdriver with SR 20 (CS 3848), easy to identify at the colored red ring, is inserted into the red T-handle with torque limitation, 6 Nm (CS 7148). The sleeve of the T-handle must be pushed forward during insertion (orange arrow in the figure).

Caution: Observe the instructions for use for the specifications for testing the torque of the instrument CS 7148.

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Securing the instrumentation

The internal locking screw is secured with the torque limiter (CS 3848 + CS 7148) with 6 Nm. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3846 + CS 7144) is used.

If the screwdriver is fully seated in the inner grub screw of the locking screw, the colored red ring then shows the correct position of the locking screw in the head of the screw.

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
Information: Use of a cross-link can be found in chapter 5.5

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Use: For reduction/restoring the sagittal or axial profile.

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Assembly of the long-arm screw and screwdriver

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The protection sleeve for longarm screws (CS 3825-03) is put over the screwdriver, open (CS 3825-01) and pushed until the sleeve clicks into place.

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By pressing the button on the retractor, see detail [2], the retractor can be moved, pushed onto the shaft wheel of the screwdriver and thus the screwdriver can be secured against unintentional release.

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Implanting the screw

The uCentum™ screw can now be inserted with the screwdriver.

Important: Hold the guide wire in position to avoid unintentional advancement. Use of a C-arm is recommended (lateral view).

To release the tight connection between the screwdriver and the screw head, the button on the retractor must be first released by pressing it and pushing it back. Then, using the shaft wheel, the screwdriver can be twisted from the screw head.

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Implanting the screw with the guide wire holder (optional)

The guidewire holder prevents the guidewire from shifting during screw insertion. When correctly used, it is an important safety feature of the system. The use of the guidewire holder is recommended particularly in the case of bicortical screw placement.

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Assembly of the guide wire holder

Important: The guidewire holder can only be used with the gearshift handle (CS 7138 or CS 7139). Use with the T-handle with ratchet (CS 7140) is not possible!

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The threaded element (CS 3822-01) is screwed into the gearshift handle as far as it will go [1]. Then the handle for the guide wire holder (CS 3822-02) is screwed into the threaded element; in doing so, the handle should not be screwed all the way into the threaded element.

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Screwing in the long-arm screw with the guide wire holder (optional)

The uCentum™ long-arm screw is placed via the guidewire in place using the screwdriver and mounted guide wire holder.

As soon as the tip of the screw comes into contact with bone, the guidewire is connected to the threaded element (CS 3822-01) by screwing in the handle (CS 3822-02) further [1].

Important: While the screw is screwed into the pedicle, the handle of the guidewire holder must be held since it should not simultaneously be rotated [2]. In this way, the guidewire is held in position while the screw is screwed in [3].

To release the tight connection between the screwdriver and the screw head, the button on the retractor must be first released by pressing it and pushing it back. Then, using the shaft wheel, the screwdriver can be twisted from the screw head.

Alignment of the screw heads

The screw heads of the screws are aligned prior to inserting the rod.


For alignment, screw head adjuster CS 3827 can be used.

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Determining rod length

The rod length can be determined using the measuring caliper (CS 5788).

Alternatively, the length and shape of the rod can be determined using the phantom rod (CS 8032-XX).

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

Inserting the rod

The uCentum™ rod is inserted using the rod holder (CS 3024). Alternatively, the rod clamp, adjustable (CS 8038) can be used for this purpose.

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

Screwing in the locking screw

Important: For all pre-fixation screws, locking screw for screws, pre-fixation CS 3801-02, must be used.

The locking screw (CS 3801-02) is screwed into the screw using the self-retaining screwdriver for pre-fixation screws (CS 3835) through the sleeve for locking screw (CS 3836) and it is gently tightened.

If necessary, the polyaxial function is completely restored upon releasing the locking screw.

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
Posture Correction

The extended screw head with threading allows for even posture correction to be performed.

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
Assembly of the counterpart

The handle for counterpart (CS 7144) must be screwed onto the counterpart, open (CS 3846).

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Clamping the polyaxial motion

The uCentum™ locking screw (CS 3801-02) is secured with the torque limiter (CS 3846 + CS 7146) with 9 Nm. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3846 + CS 7144) is used.

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If the screwdriver is fully seated in the locking screw, the colored black ring then shows the correct position of the locking screw in the head of the screw.

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
B) Alternative expansion and compression instruments (option for variable distances from 6 to 160 mm)

Selection of blades for expansion and compression

Five different blades are available. The blades can be freely combined. An expansion or compression width from 6 to 160 mm can be achieved.

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

Inserting the handle in the distraction or compression forceps

For expansion or compression, the handles (CS 3032-00) are inserted in the distraction forceps (CS 3032-01) or compression forceps (CS 3032-02). While inserting or removing the handles, press the button on the instrument (orange arrow in figure).

Push the button on the instrument to also insert or remove the blade (CS 3032-xx).

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
Parallel Distraction

For distraction, the single-part (CS 3852-02) or assembled distraction forceps (CS 3032-xx) are used. Expansion is performed gradually in the case of the single-part forceps. By pulling back the toothed rod, the forceps can be once again loosened and removed.

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
Parallel Compression

For the compression, the single-part or assembled compression forceps (CS 3852-01 or CS 3032-02 + CS 3032-00 with the corresponding blades CS 3032-05, -07, -08, -09, -10) are used. Compression is performed gradually in the case of the single-part forceps. By pulling back the toothed rod, the forceps can be once again loosened and removed.

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

Securing the instrumentation

The internal locking screw is secured with the torque limiter (CS 3848 + CS 7148) with 6 Nm. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3846 + CS 7144) is used. If the screwdriver is fully seated in the inner grub screw of the locking screw, the colored red ring then shows the correct position of the locking screw in the head of the screw.

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
Separating the long-arm blades


The blades of the long-arm screw are separated using the longarm cutting forceps (CS 8045).

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

Determining the correct implant size

The offset of both rods can be determined using the measuring caliper (CS 5788).

Important: The offset of both rods is determined on the external side of the rods [1].

The size images of the crosslinks facilitate the selection of the correct implant.

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
Implantation of the cross-link

The selected cross-link is clipped on to the two instrumentation rods using gentle pressure.

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


Securing the assembly

The locking screws are secured with the torque limiter (CS 3848 + CS 7148) with 6 Nm.

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
It should be ensured that first the two external locking screws of the cross-link are tightened and that the cross-link is thus tightly connected with the rods. Then the centrally located locking screw with the torque limiter is tightened and thus the entire cross-link is secured.

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Important: It should be ensured that the torque limiter is kept at less than a 90° angle to the rod! [1]

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





The reduction instrument is used to perform posture corrections [1].

Reduction

The inner sleeve (CS 3838-02) is positioned on the screw head [2], [3].

Then the outer sleeve (CS 3838-01) must be pushed onto the inner sleeve [4].

To facilitate screwing the outer sleeve in and out, the star grip handle (CS 3839) [5] or the socket wrench (CS 3840) [6] can be used in connection with a handle with an exchange coupling (CS 7138, CS 7139 or CS 7140).

The rod is pushed into the screw head by rotating in a clockwise direction on the hexagon of the outer sleeve. The rod is correctly positioned in the screw head when both blue instrument markings are lined up [7].

Securing the instrumentation

The locking screw is screwed into the screw head with a fully lowered outer sleeve using the screwdriver (CS 3835 or CS 3851-01) [8].


After the outer sleeve (CS 3838-01) has been removed, the inner sleeve can be removed by briefly rotating counterclockwise [9].

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In the uCentum™ system all screws except for screws with Ø 4.5 mm are provided with perforation for cement augmentation.

Mounting the adapter

Cement augmentation via the uCentum™ screw requires that a guide wire be used. It is used to facilitate assembly and ensures the patency of the screw cannulation.

The adapter (CS 3850-01) is placed on the screw via the guide wire and secured in the screw head. Then the guide wire can be removed. All adapters should be secured to the screws before starting the augmentation.

Applying the cement

Mix the cement according to the instructions and draw it into the syringes. The syringe is secured to the adapter via the Luer lock end piece. Ideally, all syringes are positioned on the screws to be augmented. Apply the cement after the waiting time has elapsed and the cement has reached the correct viscosity. The cement may need to be applied quickly depending on how fast it dries.

Important: Use a C-arm. The adapter is intended only for one-time use.

To prevent cement back flow, leave the adapters on the screws until the cement has dried. Then the adapter is released using the handle for injection adapter (CS 3850-03) and removed.

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Removal of the implant

The uCentum™ spinal system is removed in the sequence opposite to the sequence of implantation steps.

First, using the locking screws, it must be identified which screw variants, with or without a pre-fixation function, are implanted. Then the T-handle with torque limitation, 9 Nm (CS 7146) is assembled with the corresponding screwdriver (CS 3851-02 for screws without; CS 3847 for screws with a pre-fixation function). Then the torque screwdriver is used in combination with the counterpart (CS 3846 and CS 7144) to remove the locking screws.

After removing the rod, the screw is removed using the screwdriver (CS 3020-01) and the T-handle with ratchet (CS 7140).

For emergency cases, the uRevision™ Set explantation system is available.

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

Layer for uCentum screws,
polyaxial, perforated,
crosslinks and locking
screws

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Layer for uCentum rods,
length 30-180 mm,
open approach

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Layer 1.1 for uCentum
instruments open and
percutaneous approach

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
Layer 1.2 for uCentum
instruments open and
percutaneous approach

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
Layer 2.1 for uCentum
instruments open approach

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
Layer 2.2 for uCentum
instruments open approach

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Layer for uCentum screws,
pre-fixation, perforated
and locking screws

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
Layer for uCentum screws,
pre-fixation, Ø 8.5 mm and
Ø 10.0 mm, perforated

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
Layer for uCentum screws,
monoaxial, perforated
and locking screws

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Layer 1 for uCentum
longarm-screws
Ø 5.5 mm and Ø 7.5 mm,
perforated

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Layer 2 for uCentum
longarm-screws
Ø 6.5 mm, perforated

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
Layer 3.1 for
distraction and compression
open approach

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Layer 3.2 for
distraction and compression
open approach

Implantate Implants		Art.-Nr. Product number
Schraube, polyaxial, Ø 4,5 mm, Länge 25 mm	Screw, polyaxial, Ø 4.5 mm, length 25 mm	CS 3802-045-025
Schraube, polyaxial, Ø 4,5 mm, Länge 30 mm	Screw, polyaxial, Ø 4.5 mm, length 30 mm	CS 3802-045-030
Schraube, polyaxial, Ø 4,5 mm, Länge 35 mm	Screw, polyaxial, Ø 4.5 mm, length 35 mm	CS 3802-045-035
Schraube, polyaxial, Ø 4,5 mm, Länge 40 mm	Screw, polyaxial, Ø 4.5 mm, length 40 mm	CS 3802-045-040
Schraube, polyaxial, Ø 4,5 mm, Länge 45 mm	Screw, polyaxial, Ø 4.5 mm, length 45 mm	CS 3802-045-045
Schraube, polyaxial, Ø 4,5 mm, Länge 50 mm	Screw, polyaxial, Ø 4.5 mm, length 50 mm	CS 3802-045-050
Schraube, polyaxial, Ø 5,5 mm, Länge 25 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 25 mm, perforated	CS 3802-055-025
Schraube, polyaxial, Ø 5,5 mm, Länge 30 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 30 mm, perforated	CS 3802-055-030
Schraube, polyaxial, Ø 5,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 35 mm, perforated	CS 3802-055-035
Schraube, polyaxial, Ø 5,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 40 mm, perforated	CS 3802-055-040
Schraube, polyaxial, Ø 5,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 45 mm, perforated	CS 3802-055-045
Schraube, polyaxial, Ø 5,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 50 mm, perforated	CS 3802-055-050
Schraube, polyaxial, Ø 5,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 55 mm, perforated	CS 3802-055-055
Schraube, polyaxial, Ø 6,5 mm, Länge 30 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 30 mm, perforated	CS 3802-065-030
Schraube, polyaxial, Ø 6,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 35 mm, perforated	CS 3802-065-035
Schraube, polyaxial, Ø 6,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 40 mm, perforated	CS 3802-065-040
Schraube, polyaxial, Ø 6,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 45 mm, perforated	CS 3802-065-045
Schraube, polyaxial, Ø 6,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 50 mm, perforated	CS 3802-065-050
Schraube, polyaxial, Ø 6,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 55 mm, perforated	CS 3802-065-055
Schraube, polyaxial, Ø 6,5 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 60 mm, perforated	CS 3802-065-060
Schraube, polyaxial, Ø 7,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 35 mm, perforated	CS 3802-075-035
Schraube, polyaxial, Ø 7,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 40 mm, perforated	CS 3802-075-040
Schraube, polyaxial, Ø 7,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 45 mm, perforated	CS 3802-075-045
Schraube, polyaxial, Ø 7,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 50 mm, perforated	CS 3802-075-050
Schraube, polyaxial, Ø 7,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 55 mm, perforated	CS 3802-075-055
Schraube, polyaxial, Ø 7,5 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 60 mm, perforated	CS 3802-075-060
Schraube, polyaxial, Ø 8,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 35 mm, perforated	CS 3802-085-035
Schraube, polyaxial, Ø 8,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 40 mm, perforated	CS 3802-085-040

Implantate Implants		Art.-Nr. Product number
Schraube, polyaxial, Ø 8,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 45 mm, perforated	CS 3802-085-045
Schraube, polyaxial, Ø 8,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 50 mm, perforated	CS 3802-085-050
Schraube, polyaxial, Ø 8,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 55 mm, perforated	CS 3802-085-055
Schraube, polyaxial, Ø 8,5 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 60 mm, perforated	CS 3802-085-060
Schraube, polyaxial, Ø 10,0 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 45 mm, perforated	CS 3802-100-045
Schraube, polyaxial, Ø 10,0 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 50 mm, perforated	CS 3802-100-050
Schraube, polyaxial, Ø 10,0 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 55 mm, perforated	CS 3802-100-055
Schraube, polyaxial, Ø 10,0 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 60 mm, perforated	CS 3802-100-060
Schraube, monoaxial, Ø 4,5 mm, Länge 25 mm	Screw, monoaxial, Ø 4.5 mm, length 25 mm	CS 3803-045-025
Schraube, monoaxial, Ø 4,5 mm, Länge 30 mm	Screw, monoaxial, Ø 4.5 mm, length 30 mm	CS 3803-045-030
Schraube, monoaxial, Ø 4,5 mm, Länge 35 mm	Screw, monoaxial, Ø 4.5 mm, length 35 mm	CS 3803-045-035
Schraube, monoaxial, Ø 4,5 mm, Länge 40 mm	Screw, monoaxial, Ø 4.5 mm, length 40 mm	CS 3803-045-040
Schraube, monoaxial, Ø 4,5 mm, Länge 45 mm	Screw, monoaxial, Ø 4.5 mm, length 45 mm	CS 3803-045-045
Schraube, monoaxial, Ø 4,5 mm, Länge 50 mm	Screw, monoaxial, Ø 4.5 mm, length 50 mm	CS 3803-045-050
Schraube, monoaxial, Ø 5,5 mm, Länge 25 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 25 mm, perforated	CS 3803-055-025
Schraube, monoaxial, Ø 5,5 mm, Länge 30 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 30 mm, perforated	CS 3803-055-030
Schraube, monoaxial, Ø 5,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 35 mm, perforated	CS 3803-055-035
Schraube, monoaxial, Ø 5,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 40 mm, perforated	CS 3803-055-040
Schraube, monoaxial, Ø 5,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 45 mm, perforated	CS 3803-055-045
Schraube, monoaxial, Ø 5,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 50 mm, perforated	CS 3803-055-050
Schraube, monoaxial, Ø 5,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 55 mm, perforated	CS 3803-055-055
Schraube, monoaxial, Ø 6,5 mm, Länge 30 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 30 mm, perforated	CS 3803-065-030
Schraube, monoaxial, Ø 6,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 35 mm, perforated	CS 3803-065-035
Schraube, monoaxial, Ø 6,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 40 mm, perforated	CS 3803-065-040
Schraube, monoaxial, Ø 6,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 45 mm, perforated	CS 3803-065-045
Schraube, monoaxial, Ø 6,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 50 mm, perforated	CS 3803-065-050
Schraube, monoaxial, Ø 6,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 55 mm, perforated	CS 3803-065-055
Schraube, monoaxial, Ø 6,5 mm, Länge 60 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 60 mm, perforated	CS 3803-065-060

Implantate Implants		Art.-Nr. Product number
Schraube, monoaxial, Ø 7,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 35 mm, perforated	CS 3803-075-035
Schraube, monoaxial, Ø 7,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 40 mm, perforated	CS 3803-075-040
Schraube, monoaxial, Ø 7,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 45 mm, perforated	CS 3803-075-045
Schraube, monoaxial, Ø 7,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 50 mm, perforated	CS 3803-075-050
Schraube, monoaxial, Ø 7,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 55 mm, perforated	CS 3803-075-055
Schraube, monoaxial, Ø 7,5 mm, Länge 60 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 60 mm, perforated	CS 3803-075-060
Schraube, monoaxial, Ø 8,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 35 mm, perforated	CS 3803-085-035
Schraube, monoaxial, Ø 8,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 40 mm, perforated	CS 3803-085-040
Schraube, monoaxial, Ø 8,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 45 mm, perforated	CS 3803-085-045
Schraube, monoaxial, Ø 8,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 50 mm, perforated	CS 3803-085-050
Schraube, monoaxial, Ø 8,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 55 mm, perforated	CS 3803-085-055
Schraube, monoaxial, Ø 8,5 mm, Länge 60 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 60 mm, perforated	CS 3803-085-060
Schraube, vorklemmbar, Ø 4,5 mm, Länge 25 mm	Screw, pre-fixation, Ø 4.5 mm, length 25 mm	CS 3804-045-025
Schraube, vorklemmbar, Ø 4,5 mm, Länge 30 mm	Screw, pre-fixation, Ø 4.5 mm, length 30 mm	CS 3804-045-030
Schraube, vorklemmbar, Ø 4,5 mm, Länge 35 mm	Screw, pre-fixation, Ø 4.5 mm, length 35 mm	CS 3804-045-035
Schraube, vorklemmbar, Ø 4,5 mm, Länge 40 mm	Screw, pre-fixation, Ø 4.5 mm, length 40 mm	CS 3804-045-040
Schraube, vorklemmbar, Ø 4,5 mm, Länge 45 mm	Screw, pre-fixation, Ø 4.5 mm, length 45 mm	CS 3804-045-045
Schraube, vorklemmbar, Ø 4,5 mm, Länge 50 mm	Screw, pre-fixation, Ø 4.5 mm, length 50 mm	CS 3804-045-050
Schraube, vorklemmbar, Ø 5,5 mm, Länge 25 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 25 mm, perforated	CS 3804-055-025
Schraube, vorklemmbar, Ø 5,5 mm, Länge 30 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 30 mm, perforated	CS 3804-055-030
Schraube, vorklemmbar, Ø 5,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 35 mm, perforated	CS 3804-055-035
Schraube, vorklemmbar, Ø 5,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 40 mm, perforated	CS 3804-055-040
Schraube, vorklemmbar, Ø 5,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 45 mm, perforated	CS 3804-055-045
Schraube, vorklemmbar, Ø 5,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 50 mm, perforated	CS 3804-055-050
Schraube, vorklemmbar, Ø 5,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 55 mm, perforated	CS 3804-055-055
Schraube, vorklemmbar, Ø 6,5 mm, Länge 30 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 30 mm, perforated	CS 3804-065-030
Schraube, vorklemmbar, Ø 6,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 35 mm, perforated	CS 3804-065-035
Schraube, vorklemmbar, Ø 6,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 40 mm, perforated	CS 3804-065-040

Implantate Implants		Art.-Nr. Product number
Schraube, vorklemmbar, Ø 6,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 45 mm, perforated	CS 3804-065-045
Schraube, vorklemmbar, Ø 6,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 50 mm, perforated	CS 3804-065-050
Schraube, vorklemmbar, Ø 6,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 55 mm, perforated	CS 3804-065-055
Schraube, vorklemmbar, Ø 6,5 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 60 mm, perforated	CS 3804-065-060
Schraube, vorklemmbar, Ø 7,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 35 mm, perforated	CS 3804-075-035
Schraube, vorklemmbar, Ø 7,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 40 mm, perforated	CS 3804-075-040
Schraube, vorklemmbar, Ø 7,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 45 mm, perforated	CS 3804-075-045
Schraube, vorklemmbar, Ø 7,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 50 mm, perforated	CS 3804-075-050
Schraube, vorklemmbar, Ø 7,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 55 mm, perforated	CS 3804-075-055
Schraube, vorklemmbar, Ø 7,5 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 60 mm, perforated	CS 3804-075-060
Schraube, vorklemmbar, Ø 8,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 35 mm, perforated	CS 3804-085-035
Schraube, vorklemmbar, Ø 8,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 40 mm, perforated	CS 3804-085-040
Schraube, vorklemmbar, Ø 8,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 45 mm, perforated	CS 3804-085-045
Schraube, vorklemmbar, Ø 8,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 50 mm, perforated	CS 3804-085-050
Schraube, vorklemmbar, Ø 8,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 55 mm, perforated	CS 3804-085-055
Schraube, vorklemmbar, Ø 8,5 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 60 mm, perforated	CS 3804-085-060
Schraube, vorklemmbar, Ø 10,0 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 45 mm, perforated	CS 3804-100-045
Schraube, vorklemmbar, Ø 10,0 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 50 mm, perforated	CS 3804-100-050
Schraube, vorklemmbar, Ø 10,0 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 55 mm, perforated	CS 3804-100-055
Schraube, vorklemmbar, Ø 10,0 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 60 mm, perforated	CS 3804-100-060
Langkopfschraube, vorklemmbar, Ø 5,5 mm, Länge 25 mm, perforiert	Longarm screw, pre-fixation, Ø 5.5 mm, length 25 mm, perforated	CS 3805-055-025
Langkopfschraube, vorklemmbar, Ø 5,5 mm, Länge 30 mm, perforiert	Longarm screw, pre-fixation, Ø 5.5 mm, length 30 mm, perforated	CS 3805-055-030
Langkopfschraube, vorklemmbar, Ø 5,5 mm, Länge 35 mm, perforiert	Longarm screw, pre-fixation, Ø 5.5 mm, length 35 mm, perforated	CS 3805-055-035
Langkopfschraube, vorklemmbar, Ø 5,5 mm, Länge 40 mm, perforiert	Longarm screw, pre-fixation, Ø 5.5 mm, length 40 mm, perforated	CS 3805-055-040
Langkopfschraube, vorklemmbar, Ø 5,5 mm, Länge 45 mm, perforiert	Longarm screw, pre-fixation, Ø 5.5 mm, length 45 mm, perforated	CS 3805-055-045
Langkopfschraube, vorklemmbar, Ø 5,5 mm, Länge 50 mm, perforiert	Longarm screw, pre-fixation, Ø 5.5 mm, length 50 mm, perforated	CS 3805-055-050
Langkopfschraube, vorklemmbar, Ø 5,5 mm, Länge 55 mm, perforiert	Longarm screw, pre-fixation, Ø 5.5 mm, length 55 mm, perforated	CS 3805-055-055
Langkopfschraube, vorklemmbar, Ø 6,5 mm, Länge 30 mm, perforiert	Longarm screw, pre-fixation, Ø 6.5 mm, length 30 mm, perforated	CS 3805-065-030

Implantate Implants		Art.-Nr. Product number
Langkopfschraube, vorklemmbar, Ø 6,5 mm, Länge 35 mm, perforiert	Longarm screw, pre-fixation, Ø 6.5 mm, length 35 mm, perforated	CS 3805-065-035
Langkopfschraube, vorklemmbar, Ø 6,5 mm, Länge 40 mm, perforiert	Longarm screw, pre-fixation, Ø 6.5 mm, length 40 mm, perforated	CS 3805-065-040
Langkopfschraube, vorklemmbar, Ø 6,5 mm, Länge 45 mm, perforiert	Longarm screw, pre-fixation, Ø 6.5 mm, length 45 mm, perforated	CS 3805-065-045
Langkopfschraube, vorklemmbar, Ø 6,5 mm, Länge 50 mm, perforiert	Longarm screw, pre-fixation, Ø 6.5 mm, length 50 mm, perforated	CS 3805-065-050
Langkopfschraube, vorklemmbar, Ø 6,5 mm, Länge 55 mm, perforiert	Longarm screw, pre-fixation, Ø 6.5 mm, length 55 mm, perforated	CS 3805-065-055
Langkopfschraube, vorklemmbar, Ø 6,5 mm, Länge 60 mm, perforiert	Longarm screw, pre-fixation, Ø 6.5 mm, length 60 mm, perforated	CS 3805-065-060
Langkopfschraube, vorklemmbar, Ø 7,5 mm, Länge 35 mm, perforiert	Longarm screw, pre-fixation, Ø 7.5 mm, length 35 mm, perforated	CS 3805-075-035
Langkopfschraube, vorklemmbar, Ø 7,5 mm, Länge 40 mm, perforiert	Longarm screw, pre-fixation, Ø 7.5 mm, length 40 mm, perforated	CS 3805-075-040
Langkopfschraube, vorklemmbar, Ø 7,5 mm, Länge 45 mm, perforiert	Longarm screw, pre-fixation, Ø 7.5 mm, length 45 mm, perforated	CS 3805-075-045
Langkopfschraube, vorklemmbar, Ø 7,5 mm, Länge 50 mm, perforiert	Longarm screw, pre-fixation, Ø 7.5 mm, length 50 mm, perforated	CS 3805-075-050
Langkopfschraube, vorklemmbar, Ø 7,5 mm, Länge 55 mm, perforiert	Longarm screw, pre-fixation, Ø 7.5 mm, length 55 mm, perforated	CS 3805-075-055
Langkopfschraube, vorklemmbar, Ø 7,5 mm, Länge 60 mm, perforiert	Longarm screw, pre-fixation, Ø 7.5 mm, length 60 mm, perforated	CS 3805-075-060
Langkopfschraube, vorklemmbar, Ø 8,5 mm, Länge 35 mm, perforiert	Longarm screw, pre-fixation, Ø 8.5 mm, length 35 mm, perforated	CS 3805-085-035
Langkopfschraube, vorklemmbar, Ø 8,5 mm, Länge 40 mm, perforiert	Longarm screw, pre-fixation, Ø 8.5 mm, length 40 mm, perforated	CS 3805-085-040
Langkopfschraube, vorklemmbar, Ø 8,5 mm, Länge 45 mm, perforiert	Longarm screw, pre-fixation, Ø 8.5 mm, length 45 mm, perforated	CS 3805-085-045
Langkopfschraube, vorklemmbar, Ø 8,5 mm, Länge 50 mm, perforiert	Longarm screw, pre-fixation, Ø 8.5 mm, length 50 mm, perforated	CS 3805-085-050
Langkopfschraube, vorklemmbar, Ø 8,5 mm, Länge 55 mm, perforiert	Longarm screw, pre-fixation, Ø 8.5 mm, length 55 mm, perforated	CS 3805-085-055
Langkopfschraube, vorklemmbar, Ø 8,5 mm, Länge 60 mm, perforiert	Longarm screw, pre-fixation, Ø 8.5 mm, length 60 mm, perforated	CS 3805-085-060
Fixierschraube	Locking screw	CS 3801-01
Fixierschraube für Schraube, vorklemmbar	Locking screw for screw, pre-fixation	CS 3801-02
Stab, Ø 6,0 mm, gerade, Länge 30 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 30 mm, standard approach	CS 3807-030
Stab, Ø 6,0 mm, gerade, Länge 35 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 35 mm, standard approach	CS 3807-035
Stab, Ø 6,0 mm, gerade, Länge 40 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 40 mm, standard approach	CS 3807-040
Stab, Ø 6,0 mm, gerade, Länge 45 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 45 mm, standard approach	CS 3807-045
Stab, Ø 6,0 mm, gerade, Länge 50 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 50 mm, standard approach	CS 3807-050
Stab, Ø 6,0 mm, gerade, Länge 55 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 55 mm, standard approach	CS 3807-055
Stab, Ø 6,0 mm, gerade, Länge 60 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 60 mm, standard approach	CS 3807-060
Stab, Ø 6,0 mm, gerade, Länge 65 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 65 mm, standard approach	CS 3807-065

Implantate Implants		Art.-Nr. Product number
Stab, Ø 6,0 mm, gerade, Länge 70 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 70 mm, standard approach	CS 3807-070
Stab, Ø 6,0 mm, gerade, Länge 75 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 75 mm, standard approach	CS 3807-075
Stab, Ø 6,0 mm, gerade, Länge 80 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 80 mm, standard approach	CS 3807-080
Stab, Ø 6,0 mm, gerade, Länge 85 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 85 mm, standard approach	CS 3807-085
Stab, Ø 6,0 mm, gerade, Länge 90 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 90 mm, standard approach	CS 3807-090
Stab, Ø 6,0 mm, gerade, Länge 95 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 95 mm, standard approach	CS 3807-095
Stab, Ø 6,0 mm, gerade, Länge 100 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 100 mm, standard approach	CS 3807-100
Stab, Ø 6,0 mm, gerade, Länge 110 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 110 mm, standard approach	CS 3807-110
Stab, Ø 6,0 mm, gerade, Länge 120 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 120 mm, standard approach	CS 3807-120
Stab, Ø 6,0 mm, gerade, Länge 132 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 132 mm, standard approach	CS 3807-132
Stab, Ø 6,0 mm, gerade, Länge 144 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 144 mm, standard approach	CS 3807-144
Stab, Ø 6,0 mm, gerade, Länge 156 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 156 mm, standard approach	CS 3807-156
Stab, Ø 6,0 mm, gerade, Länge 168 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 168 mm, standard approach	CS 3807-168
Stab, Ø 6,0 mm, gerade, Länge 180 mm, Standardzugang	Rod, Ø 6.0 mm, straight, length 180 mm, standard approach	CS 3807-180
Querstabi variabel, XXS, Länge 35 - 43 mm	Crosslink variable, XXS, length 35 - 43 mm	CS 3810-00
Querstabi variabel, XS, Länge 41 - 56 mm	Crosslink variable, XS, length 41 - 56 mm	CS 3810-01
Querstabi variabel, S, Länge 53 - 68 mm	Crosslink variable, S, length 53 - 68 mm	CS 3810-02
Querstabi variabel, M, Länge 66 - 80 mm	crosslink variable, M, length 66 - 80 mm	CS 3810-03
Querstabi variabel, L, Länge 78 - 94 mm	Crosslink variable, L, length 78 - 94 mm	CS 3810-04
Querstabi variabel, XL, Länge 92 - 107 mm	Crosslink variable, XL, length 92 - 107 mm	CS 3810-05

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