uCentum™

comprehensive posterior system

Perkutan
Percutaneous

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+ Universal pedicle screw system
for all applications

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

Spinal systems by ulrich medical® stand for quality “Made in Germany.”

They are the successful result of systematic development activities and many years of experience in medical technology. For more than 100 years, we have done our best every day for our customers and for patient health with our proprietary and innovative products.

Leading-edge technology and competence from one source

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The small pictures at the bottom of the page show the step-by-step application of the instruments that are used as per the surgical steps on the double page. Pictures with instruments that had been used before are blue-colored.

Art. Nr./Art. No.

Art. Nr./Art. No.

This surgical technique describes the implant, the instruments, and the steps involved in the percutaneous application of the uCentum™ system. It is not sufficient as the sole basis for a successful application of the system. It is recommended to master the surgical technique with an experienced surgeon. Please observe the instructions for use, the assembly and disassembly instructions and the instrument processing instructions at www.ulrichmedical.com.

Each uCentum™ screw is cannulated and thus permits secure and precise placement via guide wires as well as minimally invasive and percutaneous application. Perforations at the tip of the screw enable cement augmentation through the already implanted screw (except for screws with Ø 4.5 mm).

The uCentum™ system is distinguished by a slim, low-profile screw head. This provides a good view for fusion monitoring.

As far as possible, the same instruments are used for polyaxial, monoaxial and pre-fixation screws. All screw variations can be freely combined. Long-arm screws and crosslinks cannot be used in percutaneous application.

The polyaxial design of the pre-fixation screws can be blocked by tightening the external locking screw with torque without, however, inhibiting any desired movement that is axial or radial to the rod. In this way, there are numerous options for reducing the spine, e.g. alignment following fractures.

The monoaxial uCentum™ screw has a higher degree of angular stability as compared to the polyaxial screws. This increased stability allows shorter instrumentation procedures and simplified reduction.

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Intended use

uCentum™ is an implant system with polyaxial, monoaxial, pre-fixation and long-arm screws, straight and preformed rods and cross stabilizers for surgical stabilization, anchorage and correction of the thoracic and lumbar spine.

Indications

Degenerative disc disease, deformity, spinal tumor, spondylolisthesis, pseudarthrosis, spinal stenosis, destruction and instability.

Contraindications

- Patients with fever or leukocytosis with non-spine associated infections
- Patients with a history of metal allergy or who tend to react to foreign bodies
- Patients whose general medical or psychological condition is unfavorable for- or could be worsened by the procedure; careful consideration is required on the part of the treating physician for these patients
- Pregnancy

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
Locking screw for screw, pre-fixation

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Rod, Ø 6.0 mm, straight
lengths 30 mm – 500 mm,
percutaneous approach

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Rod, Ø 6.0 mm, curved
lengths 30 mm – 80 mm,
percutaneous approach

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Pedicle screws

The various screw types as well as their diameters in the uCentum™ implant system are differentiated by color-coding of the screws. The screw types are differentiated by the position of the color-coding on the screw, while the diameter of the lower parts of the screw is marked with the selected color.

Caution: It should be ensured that the correct locking screw (CS 3801-01 or CS 3801-02) is used.

All pedicle screws are cannulated. They need to be implanted with the aid of a guide wire.

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CS 3801-01

Application options:

- Open and percutaneous application
- Variable-angle reduction along the rod

3

Ø 4,5 mm
Ø 4.5 mm

Ø 5,5 mm
Ø 5.5 mm

Ø 6,5 mm
Ø 6.5 mm

Ø 7,5 mm
Ø 7.5 mm

Ø 8,5 mm
Ø 8.5 mm

Ø 10 mm
Ø 10 mm


Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Ansicht im Sieb:
View in the tray:

Schraube, monoaxial | Screw, monoaxial

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





CS 3804-xxx-xxx

Position of the color-coding:
Clamping plate

Associated locking screw:
CS 3801-02

Application options:

- Open and percutaneous application
- Stable-angle reduction along the rod

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Awl, cannulated, for screws
In four diameters

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Awl, cannulated,
for screws Ø 10.0 mm,
coupling A

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Tap, cannulated,
for screws Ø 5.5 – 10.0 mm,
coupling A

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Dilator S and gauge for
screw length

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Threaded element for
guidewire holder

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Screwdriver, percutaneous,
cannulated, hex 3.5 mm

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
Screwdriver for locking
screw, self-retaining

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Rod pusher with spring,
percutaneous

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Screwdriver for locking
screw, pre-fixation, self-
retaining

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Rod pusher with spring,
percutaneous, for screws,
pre-fixation

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Sleeve for rod pusher,
percutaneous

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Socket wrench, hex 17 mm,
with coupling A

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T-handle with torque
limitation, 9 Nm
coupling V

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Screwdriver for locking
screw, coupling V

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Screwdriver for locking screw, pre-fixation, coupling V

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T-handle with torque limitation, 6 Nm, coupling H

CS 3848

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Screwdriver with SR 20, coupling H

CS 3850-02

Injektionsadapter, perkutan (Totvolumen 0,6 ml)

Adapter for injection, percutaneous (dead volume 0.6 ml)

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Handling the instrument set and the introduction of the implants are shown using a model of the spine.

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Skin incisions

A longitudinal incision of 1.5 – 2 cm is made through the skin and fascia.

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Assembling the cannulated awls

The gearshift handle (CS 7138) is positioned on the cannulated awl for screws (CS 3018-xx, CS 3823-100) and the trocar wire (CS 3821-01) is screwed in. For the various screw diameters 4.5/5.5 mm, 6.5 mm, 7.5 mm, 8.5 mm and 10 mm, there are appropriate cannulated awls available (CS 3018-04, -05, -06, -07, CS 3823-100).

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Preparation of the pedicle

The cortical bone of the pedicle is opened at the entry point using the preassembled cannulated awl for screws. The screw hole is prepared further whereby the awl is driven further through the pedicle into the vertebral body using constant rotation until it reaches the later intended position of the tip of the screw [1]. **Lateral bending of the awl should be avoided.** The trocar wire is subsequently removed [2] and the awl remains in the pedicle.

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

Inserting the guide wire

The guide wire (CS 2624) is introduced through the cannulated awl for screws [1]. The awl is then carefully removed [2] while the guide wire remains in the pedicle [3].

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Determining screw length

The gauge for screw length (CS 3820-01) is slid over the guide wire up to the bone [4]. Using the laser marking on the guide wire, the length of the screw can be read on the gauge.

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Dilatation

The additional dilators are now introduced via the guide wire in order to create a working channel for the introduction of the screws [5]. Dilator M (CS 3820-02) and dilator L (CS 3820-03) are positioned in succession. For easier placement of dilator L, the star grip handle (CS 3839) is recommended [6].

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Creating the Working Channel

Dilators S and M (CS 3820-01 and CS 3820-02) are now removed once again. Dilator L (CS 3820-03) remains in situ, also during C-arm use due to its radiolucent features.

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Tapping (optional)

As support for the pedicle preparation in the case of sclerotic cancellous vertebral bone, the tap (CS 3824-xxx) can be used. The length of the tapping instrument area allows preparation of the screw hole up to the opposing cortical bone.

Important: The guide wire must be held in position so that it is not pushed forward.

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Important: Tapping using the guide wire is not possible for screws with a diameter of 4.5 mm (CS 3824-045-xxx).

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
■ Instrumentation with polyaxial and monoaxial screws

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
Assembly of the Guide Sleeve

The stamp for guide sleeve (CS 3830-02) is placed in the outer sleeve for guide sleeve (CS 3830-01) [1] and screwed in using the assembling instrument for guide sleeve (CS 3830-04) [2]. The stamp is properly positioned when the ring on the stamp is directly on the marking on the outer sleeve [3].

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Assembly of the Screw and Guide Sleeve

The uCentum screw is introduced into the guide sleeve, aligned, as far as it will go [1]. By pressing together of both shanks and rotating the guide sleeve 90° [2], the screw audibly clicks into place in the guide sleeve [3]. In doing so, it should be ensured that the adapters on the guide sleeve fit into the grooves of the screw [4]. In the process, the screw must be kept in position. Images [5] to [8] alternatively show the use of the assembling support for guide sleeve (CS 3830-05) to facilitate the insertion of the adapters of the guide sleeve into the grooves of the screw.

First the assembly support is introduced between both shanks of the guide sleeve. During following introduction of the screw into the guide sleeve, as far as it will go, and rotating 90°, pressing together of both shanks can be dropped.

The screw is secured [7] by tightening the stamp (CS 3830-02) with the assembling instrument (CS 3830-04) [9]. The markings then lie as shown in image [10].

Note: Remove the assembling support prior to implantation, and check the tight connection between the screw and guide sleeve.

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Assembly of the Screw and Screwdriver

The screwdriver (CS 3826) is inserted into the guide sleeve and screwed in fully. The T-handle (CS 7140) or gearshift handle (CS 7138 / 7139) may be used as a handle.

Note: In the case of polyaxial screws (CS 3802-xxx-xxx) and prefixation screws (CS 3804-xxx-xxx), the screw head does not rotate with the guide sleeve when the screw is screwed in, but during use of monoaxial screws.

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Implanting the Screw without Guide Wire Holder

The uCentum screw is inserted via the guide wire in place using the screwdriver (CS 3826), including guide sleeve.

Important: The guide wire must be held in position so that it is not pushed forward when the screw is screwed in. Lateral C-arm control.

To loosen the tight connection, the screwdriver must be released from the guide sleeve.

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Instrumentation with polyaxial and monoaxial screws

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Implanting the Screw with the Guide Wire Holder

The guide wire holder (CS 3822-XX) prevents the guide wire from shifting forward during screw insertion. When correctly used, it is an important safety feature of the system. The use of the guide wire holder is recommended particularly in the case of bicortical screw placement.

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Assembly of the Guide Wire Holder

Important: The guide wire holder can only be used with the gear-shift handle (CS 7138 or CS 7139). Use with the T-handle with ratchet (CS 7140) is not possible!

The threaded element (CS 3822-01) is screwed into one of the two gear shift handles as far as it will go [1]. Then the handle for the guide wire holder (CS 3822-02) is screwed into the threaded element [2]; in doing so, the handle should not be screwed all the way into the threaded element.

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
Screwing in the Screw with the Guide Wire Holder

The uCentum screw is inserted via the guide wire in place using the screwdriver and mounted guide wire holder.

As soon as the tip of the screw comes into contact with bone, the guide wire is secured to the threaded element (CS 3822-01) by screwing in the handle (CS 3822-02) further [1].

Important: While the screw is screwed into the pedicle, the guide wire holder should not simultaneously be rotated [2]. In this way, the guide wire is held in position while the screw is screwed in [3].

In case additional screws are implanted, it should be ensured that the sleeves are parallel.

To loosen the tight connection, the screwdriver must be rotated out of the guide sleeve. The screwdriver, guide wire and dilator L (CS 3820-03) are removed; the guide sleeve remains in situ.

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Determining Rod Length

The length and shape of the rod can be determined using the phantom rod (CS 8032-XX). For this purpose, the guide sleeves must be aligned in parallel.

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

Assembly of the Rod using the Rod Inserter, S-shaped

The rod for the percutaneous approach (CS 3808-xxx and CS 3809-xxx) is inserted into the intended recess of the rod inserter, S-shaped (CS 3832) [1]. The correct alignment of the rod to the rod inserter is shown by the marking and position of the recesses on the rod [2]. A tight connection is ensured with the aid of the screwdriver (CS 3831) [3].

As needed, the rod can be brought into the suitable shape using the rod bending forceps according to Lepine (UT 1639-29) [4].

Important: Prior to implantation, check the tight connection between the rod and rod inserter.

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
Inserting the rod using the Rod Inserter, S-shaped

The uCentum™ rod is guided from a cranial or caudal direction using the rod inserter, S-shaped (CS 3832), into the opening of the first guide sleeve [1]. By gently rotating the rod inserter in the frontal plane, it can be determined whether the rod is inside the other guide sleeves [2]. Lateral C-arm control.

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
Assembly of the Rod using the Rod inserter, L-shaped

The rod for the percutaneous approach (CS 3808-xxx and CS 3809-xxx) is inserted into the intended recess of the rod inserter, L-shaped (CS 3833) [1]. The correct alignment of the rod to the rod inserter is shown by the marking and position of the recesses on the rod [2]. A tight connection is ensured with the aid of the screwdriver (CS 3831) [3].

As needed, the rod can be brought into the suitable shape using the rod bending forceps according to Lepine (UT 1639-29) [4].

Important: Prior to implantation, check the tight connection between the rod and rod inserter.

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

Inserting the Rod using the Rod Inserter, L-shaped

The uCentum™ rod is guided using the rod inserter, L-shaped (CS 3833), into the opening of the first guide sleeve [1]. By gently rotating the rod inserter in the frontal plane, it can be determined whether the rod is inside the other guide sleeves [2]. Lateral C-arm control.

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
Assembly of the Rod Inserter, curved (for short and long instrumentations)

The handle (CS 3834-03) is placed in the desired position on the tube of the rod inserter (CS 3834-01) and inserted as far as it will go [1]. Then the stamp (CS 3834-02) is introduced into the tube [2].

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
Assembly of the Rod using the Rod Inserter

The rod (CS 3808-xxx or CS 3809-xxx) is snapped into place, aligned, in the tip of the rod inserter (CS 3834) [3]. The correct alignment of the rod to the rod inserter is shown by the marking and position of the recesses on the rod. By tightening the set screw (CS 3834-02) by hand, the rod is firmly connected to the rod inserter [4].

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
As needed, the rod can be appropriately shaped using the rod bending forceps according to Lepine (UT 1639-29) [5].

Important: Prior to implantation, check the tight connection between the rod and rod inserter.

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
Inserting the Rod using the Rod Inserter

The uCentum rod is bent using the rod inserter (CS 3834) and guided via an additional incision into the opening of the first guide sleeve [1]. By gently rotating the rod inserter in the frontal plane, it can be determined whether the rod is inside the other guide sleeves. Lateral C-arm control.

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Checking Rod Position

To check whether the rod is within the head of the screw and whether there is sufficient projection, the dilator S (CS 3820-01) is introduced in the guide sleeve. The laser marking on the dilator is visible when the rod is correctly positioned.

Important: The rod should project at least 5 mm. The rod inserter should not be inside the first guide sleeve. Lateral C-arm control.

Please observe chapter 5.5. if reduction should become necessary.

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Screwing in the Locking Screw

The uCentum™ locking screw (CS 3801-01) is screwed into the screw through the guide sleeve using the self-retaining screwdriver (CS 3851-01) and gently secured. The rod and also the polyaxial head are fixed.

If necessary, the polyaxial function is completely restored upon releasing the locking screw.

Important: In the case of a reduction (see section 5.5.), the locking screw is not screwed in until the rod has been completely lowered into the screw head.

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Assembly of the Torque Limiter

The screwdriver for locking screw (CS 3851-02) is inserted into the T-handle with torque limitation (CS 7146). The sleeve of the T-handle must be pulled back during insertion (orange arrow in the figure).

Caution: Observe the instructions for use for the specifications for testing the torque of the instruments CS 7146 and CS 7148.

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Assembly of the Counterpart

The L handle, M8 (CS 7144) must be screwed onto the counterpart, percutaneous (CS 3849).

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Securing the Screw Head

The uCentum™ locking screw is secured with the torque limiter (CS 7146 + CS 3851-02) with 9 Nm. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3849 + CS 7144) is used. If the screwdriver is fully seated on the outer ring of the locking screw, the colored black ring will show the correct position of the locking screw in the screw head.

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
Releasing the Connection between the Rod and Rod Insertor

To release the rod from the rod inserter, S- and L-shaped, the screwdriver for rod inserter SR 20 (CS 3831) must be brought down laterally on the guide sleeve to the rod inserter. The fixation is released by introducing the screwdriver and turning it counter-clockwise. By rotating the set screw of the rod inserter, curved counterclockwise, the rod is disconnected. Then the rod inserter can be removed by carefully tilting and disengaging it.

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
Removal of the Guide Sleeve

After completing the assembly, the stamp (CS 3830-02) with the assembling instrument for guide sleeve (CS 3830-04) is released by one rotation counter clockwise [1]. By subsequently rotating the guide sleeve (CS 3830-01 and -02) by 90°, the connection to the screw is released with an audible click [2]. The guide sleeve can be removed.

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
■ Instrumentation with pre-fixation screws

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Screwing in the Locking Screw for Pre-Fixation Screws

The uCentum™ locking screw (CS 3801-02) is screwed into the screw through the guide sleeve using the self-retaining screwdriver (CS 3835) and gently secured.

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The final fixation takes place using the torque limiter with 9 Nm.

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
Assembly of the Torque Limiter

The screwdriver for locking screw, pre-fixation (CS 3847) is inserted into the T-handle with torque limitation (CS 7146). The sleeve of the T-handle must be pulled back during insertion (orange arrow in the figure).

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
**Securing the
Screw Head**

The uCentum™ locking screw (CS 3801-02) is secured with the torque limiter (CS 7146 + CS 3847) with 9 Nm. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3849 + CS 7144) is used. If the screwdriver is fully seated on the outer ring of the locking screw, the colored black ring will show the correct position of the locking screw in the screw head.

Important: Only the polyaxial head is secured in this way. The screw head remains movable along the rod and radially to the rod.

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
Assembly of the Torque Limiter

The screwdriver with SR 20 (CS 3848), easily recognizable by the colored red ring, is inserted into the red T-handle with torque limitation, 6 Nm (CS 7148). The sleeve of the T-handle must be pushed forward during insertion (orange arrow in the figure).

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Securing the Assembly

The uCentum™ locking screw (CS 3801-02) is secured with the T-handle with torque limitation, 6 Nm (CS 7148 + CS 3848) with 6 Nm. To relieve the instrumented section of the spinal column and for secure fixation, the assembled counterpart (CS 3849 + CS 7144) is used. If the screwdriver is fully seated in the inner grub screw of the locking screw, the colored red ring will show the correct position of the locking screw in the screw head.

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
Releasing the Connection between the Rod and Rod Inserter

To release the rod from the rod inserter, S- and L-shaped, the screwdriver for rod inserter SR 20 (CS 3831) must be brought down laterally on the guide sleeve to the rod inserter. The fixation is released by introducing the screwdriver and turning it counter-clockwise. By rotating the set screw of the rod inserter, curved counterclockwise, the rod is disconnected. Then the rod inserter can be removed by carefully tilting and dis-engaging it.

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
Reduction with sleeve and nuts

If the rod is resting above the screw head or if posture correction is to be performed, the sleeve (CS 3841-01) is pushed down over the guide sleeve to the rod [1]. Screwing in the nut (CS 3841-03) in the guide sleeve lowers the rod in the screw head [2]. The reduction channel is indicated by lines marked on the sleeve (at 5-mm intervals). When the rod is fully inserted in the screw head, the orange markings line up.



To facilitate operating the nuts, the star grip handle, hex 17 mm (CS 3839) or the socket wrench, hex 17 mm (CS 3840) can be used in connection with a handle with an exchange coupling (CS 7138 or CS 7140).

After the rod is completely positioned in the screw head the proper locking screw can be inserted with the screwdriver.

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**Reduction using Rod
Pusher with Spring,
percutaneous**

Depending on the screw type, the corresponding locking screw (CS 3801-01 /-02) is placed on the rod pusher with spring, percutaneous (CS 3829 or CS 3829-01) [1].

The rod pusher (CS 3829) is screwed in as far as it will go into the guide sleeve using a locking screw, in this way, the rod is pressed into the head of the screw [2].

To facilitate manipulating the nuts, the star grip handle, hex 17 mm (CS 3839) or the socket wrench, hex 17 mm (CS 3840) can be used in connection with a handle with an exchange coupling (CS 7138, CS 7139 or CS 7140). Then the locking screw can be screwed into the screw head using screwdriver for rod inserter, SR 20 (CS 3831) [3]. The rest of the procedure for final fixation is described in chapter 5.2. (locking screw CS 3801-01 on p. 31) and chapter 5.3 (locking screw for screw, pre-fixation CS 3801-02 on p. 46 and 47).

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
Cement augmentation

Screws with lateral perforations can be augmented with cement as needed.

Mounting the Adapter

Cement augmentation via the uCentum™ screw requires that a guide wire be used. It is used to facilitate assembly and ensures the patency of the screw cannulation.

The adapter for injection, percutaneous (CS 3850-02) is placed on the screw via the guide wire through the guide sleeve and secured in the screw head. Then the guide wire can be removed. All adapters should be secured to the screws before starting the augmentation.

Applying the cement

Mix the cement according to the instructions and draw it into the syringes. The syringe is secured to an adapter via the Luer lock end piece. Ideally, all syringes are positioned on the screws to be augmented. Apply the cement after the waiting time has elapsed and the cement has reached the correct viscosity. The cement may need to be applied quickly depending on how fast it dries.

Important: Use a C-arm. The adapter is intended only for one-time use.

To prevent cement back flow, leave the adapters on the screws until the cement has dried. Then the adapter is released using the handle for injection adapters CS 3850-03 and removed.

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
Layer for uCentum screws,
polyaxial, perforated,
crosslinks and locking
screws

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Layer for uCentum rods,
length 30-180 mm,
open approach

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Layer 1.1 for uCentum
instruments open and
percutaneous approach

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
Layer 1.2 for uCentum
instruments open and
percutaneous approach

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
Layer for uCentum rods, curved, length 30-80 mm, percutaneous approach

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
Layer for uCentum rods, straight, length 30-80 mm, percutaneous approach

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
Layer for uCentum rods, straight, length 85-400 mm, percutaneous approach

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Layer P1.1 for uCentum instruments, percutaneous approach

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
Layer P1.2 for uCentum instruments, percutaneous approach

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
Layer for uCentum screws,
pre-fixation, perforated
and locking screws

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Layer for uCentum screws,
pre-fixation, Ø 8.5 mm and
10.0 mm, perforated

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Layer for uCentum screws,
monoaxial, perforated
and locking screws

Implantate Implants		Art.-Nr. Product number
Schraube, polyaxial, Ø 4,5 mm, Länge 25 mm	Screw, polyaxial, Ø 4.5 mm, length 25 mm	CS 3802-045-025
Schraube, polyaxial, Ø 4,5 mm, Länge 30 mm	Screw, polyaxial, Ø 4.5 mm, length 30 mm	CS 3802-045-030
Schraube, polyaxial, Ø 4,5 mm, Länge 35 mm	Screw, polyaxial, Ø 4.5 mm, length 35 mm	CS 3802-045-035
Schraube, polyaxial, Ø 4,5 mm, Länge 40 mm	Screw, polyaxial, Ø 4.5 mm, length 40 mm	CS 3802-045-040
Schraube, polyaxial, Ø 4,5 mm, Länge 45 mm	Screw, polyaxial, Ø 4.5 mm, length 45 mm	CS 3802-045-045
Schraube, polyaxial, Ø 4,5 mm, Länge 50 mm	Screw, polyaxial, Ø 4.5 mm, length 50 mm	CS 3802-045-050
Schraube, polyaxial, Ø 5,5 mm, Länge 25 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 25 mm, perforated	CS 3802-055-025
Schraube, polyaxial, Ø 5,5 mm, Länge 30 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 30 mm, perforated	CS 3802-055-030
Schraube, polyaxial, Ø 5,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 35 mm, perforated	CS 3802-055-035
Schraube, polyaxial, Ø 5,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 40 mm, perforated	CS 3802-055-040
Schraube, polyaxial, Ø 5,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 45 mm, perforated	CS 3802-055-045
Schraube, polyaxial, Ø 5,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 50 mm, perforated	CS 3802-055-050
Schraube, polyaxial, Ø 5,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 5.5 mm, length 55 mm, perforated	CS 3802-055-055
Schraube, polyaxial, Ø 6,5 mm, Länge 30 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 30 mm, perforated	CS 3802-065-030
Schraube, polyaxial, Ø 6,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 35 mm, perforated	CS 3802-065-035
Schraube, polyaxial, Ø 6,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 40 mm, perforated	CS 3802-065-040
Schraube, polyaxial, Ø 6,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 45 mm, perforated	CS 3802-065-045
Schraube, polyaxial, Ø 6,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 50 mm, perforated	CS 3802-065-050
Schraube, polyaxial, Ø 6,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 55 mm, perforated	CS 3802-065-055
Schraube, polyaxial, Ø 6,5 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 6.5 mm, length 60 mm, perforated	CS 3802-065-060
Schraube, polyaxial, Ø 7,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 35 mm, perforated	CS 3802-075-035
Schraube, polyaxial, Ø 7,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 40 mm, perforated	CS 3802-075-040
Schraube, polyaxial, Ø 7,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 45 mm, perforated	CS 3802-075-045
Schraube, polyaxial, Ø 7,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 50 mm, perforated	CS 3802-075-050
Schraube, polyaxial, Ø 7,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 55 mm, perforated	CS 3802-075-055
Schraube, polyaxial, Ø 7,5 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 7.5 mm, length 60 mm, perforated	CS 3802-075-060
Schraube, polyaxial, Ø 8,5 mm, Länge 35 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 35 mm, perforated	CS 3802-085-035
Schraube, polyaxial, Ø 8,5 mm, Länge 40 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 40 mm, perforated	CS 3802-085-040
Schraube, polyaxial, Ø 8,5 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 45 mm, perforated	CS 3802-085-045

Implantate Implants		Art.-Nr. Product number
Schraube, polyaxial, Ø 8,5 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 50 mm, perforated	CS 3802-085-050
Schraube, polyaxial, Ø 8,5 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 55 mm, perforated	CS 3802-085-055
Schraube, polyaxial, Ø 8,5 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 8.5 mm, length 60 mm, perforated	CS 3802-085-060
Schraube, polyaxial, Ø 10,0 mm, Länge 45 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 45 mm, perforated	CS 3802-100-045
Schraube, polyaxial, Ø 10,0 mm, Länge 50 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 50 mm, perforated	CS 3802-100-050
Schraube, polyaxial, Ø 10,0 mm, Länge 55 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 55 mm, perforated	CS 3802-100-055
Schraube, polyaxial, Ø 10,0 mm, Länge 60 mm, perforiert	Screw, polyaxial, Ø 10,0 mm, length 60 mm, perforated	CS 3802-100-060
Schraube, monoaxial, Ø 4,5 mm, Länge 25 mm	Screw, monoaxial, Ø 4.5 mm, length 25 mm	CS 3803-045-025
Schraube, monoaxial, Ø 4,5 mm, Länge 30 mm	Screw, monoaxial, Ø 4.5 mm, length 30 mm	CS 3803-045-030
Schraube, monoaxial, Ø 4,5 mm, Länge 35 mm	Screw, monoaxial, Ø 4.5 mm, length 35 mm	CS 3803-045-035
Schraube, monoaxial, Ø 4,5 mm, Länge 40 mm	Screw, monoaxial, Ø 4.5 mm, length 40 mm	CS 3803-045-040
Schraube, monoaxial, Ø 4,5 mm, Länge 45 mm	Screw, monoaxial, Ø 4.5 mm, length 45 mm	CS 3803-045-045
Schraube, monoaxial, Ø 4,5 mm, Länge 50 mm	Screw, monoaxial, Ø 4.5 mm, length 50 mm	CS 3803-045-050
Schraube, monoaxial, Ø 5,5 mm, Länge 25 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 25 mm, perforated	CS 3803-055-025
Schraube, monoaxial, Ø 5,5 mm, Länge 30 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 30 mm, perforated	CS 3803-055-030
Schraube, monoaxial, Ø 5,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 35 mm, perforated	CS 3803-055-035
Schraube, monoaxial, Ø 5,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 40 mm, perforated	CS 3803-055-040
Schraube, monoaxial, Ø 5,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 45 mm, perforated	CS 3803-055-045
Schraube, monoaxial, Ø 5,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 50 mm, perforated	CS 3803-055-050
Schraube, monoaxial, Ø 5,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 5.5 mm, length 55 mm, perforated	CS 3803-055-055
Schraube, monoaxial, Ø 6,5 mm, Länge 30 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 30 mm, perforated	CS 3803-065-030
Schraube, monoaxial, Ø 6,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 35 mm, perforated	CS 3803-065-035
Schraube, monoaxial, Ø 6,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 40 mm, perforated	CS 3803-065-040
Schraube, monoaxial, Ø 6,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 45 mm, perforated	CS 3803-065-045
Schraube, monoaxial, Ø 6,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 50 mm, perforated	CS 3803-065-050
Schraube, monoaxial, Ø 6,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 55 mm, perforated	CS 3803-065-055
Schraube, monoaxial, Ø 6,5 mm, Länge 60 mm, perforiert	Screw, monoaxial, Ø 6.5 mm, length 60 mm, perforated	CS 3803-065-060
Schraube, monoaxial, Ø 7,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 35 mm, perforated	CS 3803-075-035
Schraube, monoaxial, Ø 7,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 40 mm, perforated	CS 3803-075-040

Implantate Implants		Art.-Nr. Product number
Schraube, monoaxial, Ø 7,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 45 mm, perforated	CS 3803-075-045
Schraube, monoaxial, Ø 7,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 50 mm, perforated	CS 3803-075-050
Schraube, monoaxial, Ø 7,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 55 mm, perforated	CS 3803-075-055
Schraube, monoaxial, Ø 7,5 mm, Länge 60 mm, perforiert	Screw, monoaxial, Ø 7.5 mm, length 60 mm, perforated	CS 3803-075-060
Schraube, monoaxial, Ø 8,5 mm, Länge 35 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 35 mm, perforated	CS 3803-085-035
Schraube, monoaxial, Ø 8,5 mm, Länge 40 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 40 mm, perforated	CS 3803-085-040
Schraube, monoaxial, Ø 8,5 mm, Länge 45 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 45 mm, perforated	CS 3803-085-045
Schraube, monoaxial, Ø 8,5 mm, Länge 50 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 50 mm, perforated	CS 3803-085-050
Schraube, monoaxial, Ø 8,5 mm, Länge 55 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 55 mm, perforated	CS 3803-085-055
Schraube, monoaxial, Ø 8,5 mm, Länge 60 mm, perforiert	Screw, monoaxial, Ø 8.5 mm, length 60 mm, perforated	CS 3803-085-060
Schraube, vorklemmbar, Ø 4,5 mm, Länge 25 mm	Screw, pre-fixation, Ø 4.5 mm, length 25 mm	CS 3804-045-025
Schraube, vorklemmbar, Ø 4,5 mm, Länge 30 mm	Screw, pre-fixation, Ø 4.5 mm, length 30 mm	CS 3804-045-030
Schraube, vorklemmbar, Ø 4,5 mm, Länge 35 mm	Screw, pre-fixation, Ø 4.5 mm, length 35 mm	CS 3804-045-035
Schraube, vorklemmbar, Ø 4,5 mm, Länge 40 mm	Screw, pre-fixation, Ø 4.5 mm, length 40 mm	CS 3804-045-040
Schraube, vorklemmbar, Ø 4,5 mm, Länge 45 mm	Screw, pre-fixation, Ø 4.5 mm, length 45 mm	CS 3804-045-045
Schraube, vorklemmbar, Ø 4,5 mm, Länge 50 mm	Screw, pre-fixation, Ø 4.5 mm, length 50 mm	CS 3804-045-050
Schraube, vorklemmbar, Ø 5,5 mm, Länge 25 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 25 mm, perforated	CS 3804-055-025
Schraube, vorklemmbar, Ø 5,5 mm, Länge 30 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 30 mm, perforated	CS 3804-055-030
Schraube, vorklemmbar, Ø 5,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 35 mm, perforated	CS 3804-055-035
Schraube, vorklemmbar, Ø 5,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 40 mm, perforated	CS 3804-055-040
Schraube, vorklemmbar, Ø 5,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 45 mm, perforated	CS 3804-055-045
Schraube, vorklemmbar, Ø 5,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 50 mm, perforated	CS 3804-055-050
Schraube, vorklemmbar, Ø 5,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 5.5 mm, length 55 mm, perforated	CS 3804-055-055
Schraube, vorklemmbar, Ø 6,5 mm, Länge 30 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 30 mm, perforated	CS 3804-065-030
Schraube, vorklemmbar, Ø 6,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 35 mm, perforated	CS 3804-065-035
Schraube, vorklemmbar, Ø 6,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 40 mm, perforated	CS 3804-065-040
Schraube, vorklemmbar, Ø 6,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 45 mm, perforated	CS 3804-065-045
Schraube, vorklemmbar, Ø 6,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 50 mm, perforated	CS 3804-065-050
Schraube, vorklemmbar, Ø 6,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 55 mm, perforated	CS 3804-065-055

Implantate Implants		Art.-Nr. Product number
Schraube, vorklemmbar, Ø 6,5 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 6.5 mm, length 60 mm, perforated	CS 3804-065-060
Schraube, vorklemmbar, Ø 7,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 35 mm, perforated	CS 3804-075-035
Schraube, vorklemmbar, Ø 7,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 40 mm, perforated	CS 3804-075-040
Schraube, vorklemmbar, Ø 7,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 45 mm, perforated	CS 3804-075-045
Schraube, vorklemmbar, Ø 7,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 50 mm, perforated	CS 3804-075-050
Schraube, vorklemmbar, Ø 7,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 55 mm, perforated	CS 3804-075-055
Schraube, vorklemmbar, Ø 7,5 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 7.5 mm, length 60 mm, perforated	CS 3804-075-060
Schraube, vorklemmbar, Ø 8,5 mm, Länge 35 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 35 mm, perforated	CS 3804-085-035
Schraube, vorklemmbar, Ø 8,5 mm, Länge 40 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 40 mm, perforated	CS 3804-085-040
Schraube, vorklemmbar, Ø 8,5 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 45 mm, perforated	CS 3804-085-045
Schraube, vorklemmbar, Ø 8,5 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 50 mm, perforated	CS 3804-085-050
Schraube, vorklemmbar, Ø 8,5 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 55 mm, perforated	CS 3804-085-055
Schraube, vorklemmbar, Ø 8,5 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 8.5 mm, length 60 mm, perforated	CS 3804-085-060
Schraube, vorklemmbar, Ø 10,0 mm, Länge 45 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 45 mm, perforated	CS 3804-100-045
Schraube, vorklemmbar, Ø 10,0 mm, Länge 50 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 50 mm, perforated	CS 3804-100-050
Schraube, vorklemmbar, Ø 10,0 mm, Länge 55 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 55 mm, perforated	CS 3804-100-055
Schraube, vorklemmbar, Ø 10,0 mm, Länge 60 mm, perforiert	Screw, pre-fixation, Ø 10.0 mm, length 60 mm, perforated	CS 3804-100-060
Fixierschraube	Locking screw	CS 3801-01
Fixierschraube für Schraube, vorklemmbar	Locking screw for screw, pre-fixation	CS 3801-02
Stab, Ø 6,0 mm, gerade, Länge 30 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 30 mm, percutaneous approach	CS 3808-030
Stab, Ø 6,0 mm, gerade, Länge 35 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 35 mm, percutaneous approach	CS 3808-035
Stab, Ø 6,0 mm, gerade, Länge 40 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 40 mm, percutaneous approach	CS 3808-040
Stab, Ø 6,0 mm, gerade, Länge 45 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 45 mm, percutaneous approach	CS 3808-045
Stab, Ø 6,0 mm, gerade, Länge 50 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 50 mm, percutaneous approach	CS 3808-050
Stab, Ø 6,0 mm, gerade, Länge 55 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 55 mm, percutaneous approach	CS 3808-055
Stab, Ø 6,0 mm, gerade, Länge 60 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 60 mm, percutaneous approach	CS 3808-060
Stab, Ø 6,0 mm, gerade, Länge 65 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 65 mm, percutaneous approach	CS 3808-065
Stab, Ø 6,0 mm, gerade, Länge 70 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 70 mm, percutaneous approach	CS 3808-070
Stab, Ø 6,0 mm, gerade, Länge 75 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 75 mm, percutaneous approach	CS 3808-075

Implantate Implants		Art.-Nr. Product number
Stab, Ø 6,0 mm, gerade, Länge 80 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 80 mm, percutaneous approach	CS 3808-080
Stab, Ø 6,0 mm, gerade, Länge 85 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 85 mm, percutaneous approach	CS 3808-085
Stab, Ø 6,0 mm, gerade, Länge 90 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 90 mm, percutaneous approach	CS 3808-090
Stab, Ø 6,0 mm, gerade, Länge 95 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 95 mm, percutaneous approach	CS 3808-095
Stab, Ø 6,0 mm, gerade, Länge 100 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 100 mm, percutaneous approach	CS 3808-100
Stab, Ø 6,0 mm, gerade, Länge 110 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 110 mm, percutaneous approach	CS 3808-110
Stab, Ø 6,0 mm, gerade, Länge 120 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 120 mm, percutaneous approach	CS 3808-120
Stab, Ø 6,0 mm, gerade, Länge 132 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 132 mm, percutaneous approach	CS 3808-132
Stab, Ø 6,0 mm, gerade, Länge 144 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 144 mm, percutaneous approach	CS 3808-144
Stab, Ø 6,0 mm, gerade, Länge 156 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 156 mm, percutaneous approach	CS 3808-156
Stab, Ø 6,0 mm, gerade, Länge 168 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 168 mm, percutaneous approach	CS 3808-168
Stab, Ø 6,0 mm, gerade, Länge 180 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 180 mm, percutaneous approach	CS 3808-180
Stab, Ø 6,0 mm, gerade, Länge 200 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 200 mm, percutaneous approach	CS 3808-200
Stab, Ø 6,0 mm, gerade, Länge 300 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 300 mm, percutaneous approach	CS 3808-300
Stab, Ø 6,0 mm, gerade, Länge 400 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 400 mm, percutaneous approach	CS 3808-400
Stab, Ø 6,0 mm, gebogen, Länge 30 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 30 mm, percutaneous approach	CS 3809-030
Stab, Ø 6,0 mm, gebogen, Länge 35 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 35 mm, percutaneous approach	CS 3809-035
Stab, Ø 6,0 mm, gebogen, Länge 40 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 40 mm, percutaneous approach	CS 3809-040
Stab, Ø 6,0 mm, gebogen, Länge 45 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 45 mm, percutaneous approach	CS 3809-045
Stab, Ø 6,0 mm, gebogen, Länge 50 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 50 mm, percutaneous approach	CS 3809-050
Stab, Ø 6,0 mm, gebogen, Länge 55 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 55 mm, percutaneous approach	CS 3809-055
Stab, Ø 6,0 mm, gebogen, Länge 60 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 60 mm, percutaneous approach	CS 3809-060
Stab, Ø 6,0 mm, gebogen, Länge 65 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 65 mm, percutaneous approach	CS 3809-065
Stab, Ø 6,0 mm, gebogen, Länge 70 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 70 mm, percutaneous approach	CS 3809-070
Stab, Ø 6,0 mm, gebogen, Länge 75 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 75 mm, percutaneous approach	CS 3809-075
Stab, Ø 6,0 mm, gebogen, Länge 80 mm, perkutaner Zugang	Rod, Ø 6.0 mm, curved, length 80 mm, percutaneous approach	CS 3809-080
Auf Anfrage optional erhältlich:		Available on request:
Stab, Ø 6,0 mm, gerade, Länge 500 mm, perkutaner Zugang	Rod, Ø 6.0 mm, straight, length 500 mm, percutaneous approach	CS 3808-500

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