

Introducing PediGuard® XS

The PediGuard® technology is now available with a Stiffer, Sharper and Miniaturized tip.

Like the original PediGuard®, the **NEW** PediGuard® XS takes the place of the standard gearshift while alerting the surgeon prior to a breach by accurately analyzing the electrical conductivity of the surrounding tissues. The surgeon is informed in real time by changes in pitch and cadence of audio feedback and a flashing LED light.




Classic
Curved
 Cannulated

PediGuard®


The **NEW** PediGuard® Curved XS is now available in a slimmer design to address small sized pedicles.

NEW

PediGuard® Curv XS
 (Ref. P1-AU451)



PediGuard® Curv
 (Ref. P1-AU450)




Classic
 Curved
Cannulated

PediGuard®

The **NEW** PediGuard® 2.5mm XS is now available in a sharper and stiffer version to enhance penetration into harder bone.

NEW

PediGuard® Tri Tip Ø2.5mm XS
 (Ref. P1-AU414)



PediGuard® Tri Tip Ø2.5mm
 (Ref. P1-AU413)



SpineGuard® S.A.
 5-7, rue de l'Amiral Courbet
 94160 Saint-Mandé - France
 Phone: +33 1 45 18 45 19
 Fax: +33 1 45 18 45 20
 Manufacturer

SpineGuard®
 Make spine surgery safer
www.spineguard.com

SpineGuard® Inc.
 1388 Sutter Street, Suite 510
 San Francisco, CA 94109 - USA
 Phone: +1 415 512 2500
 Fax: +1 415 512 8004
 Distributor

Classic
Curved
Cannulated

PediGuard®



XS

Stiffer.
Sharper.
Miniaturized



SpineGuard®
Make spine surgery safer

www.spineguard.com

CE 0123

Caution: See package insert for labeling limitations, intended uses, relevant warnings, precautions, side effects and contraindications. Federal (USA) law restricts the sale and use of this device to a prescription of a physician

Ref. GP PXS113E

NON CONTRACTUAL DOCUMENT