

ДЕНТАЛЬНАЯ ИМПЛАНТОЛОГИЯ

СПИНАЛЬНАЯ ХИРУРГИЯ

ЧЕРЕПНО-ЧЕЛЮСТНО-ЛИЦЕВАЯ
ХИРУРГИЯ

ДЕНТАЛЬНАЯ ИМПЛАНТОЛОГИЯ

СПИНАЛЬНАЯ ХИРУРГИЯ

ЧЕРЕПНО-ЧЕЛЮСТНО-ЛИЦЕВАЯ
ХИРУРГИЯ

ДЕНТАЛЬНАЯ ИМПЛАНТОЛОГИЯ

СПИНАЛЬНАЯ ХИРУРГИЯ

ДЕНТАЛЬНАЯ ИМПЛАНТОЛОГИЯ

**ЧЕРЕПНО-ЧЕЛЮСТНО-ЛИЦЕВАЯ
ХИРУРГИЯ**

СПИНАЛЬНАЯ ХИРУРГИЯ

СПИНАЛЬНАЯ ХИРУРГИЯ

ЧЕРЕПНО-ЧЕЛЮСТНО-ЛИЦЕВАЯ
ХИРУРГИЯ

ДЕНТАЛЬНАЯ ИМПЛАНТОЛОГИЯ

ЧЕРЕПНО-ЧЕЛЮСТНО-ЛИЦЕВАЯ
ХИРУРГИЯ

СПИНАЛЬНАЯ ХИРУРГИЯ

ДЕНТАЛЬНАЯ ИМПЛАНТОЛОГИЯ



Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 9001:2008

This is to certify that

Conmet, LLC
24/1 Onezhskaya str.
Moscow
125413
Russian Fed.

Holds Certificate No: **PM 809974**

and operates a Quality Management System which complies with the requirements of ISO 9001:2008 for the scope:

Design, development and production of fabricate metal products and fine metal mechanical, metal medical instruments and implants.

For and on behalf of BSI:

Gary Fenton, Global Assurance Director

Originally registered: 05/10/2008

Latest Issue: 05/07/2012

Expiry Date: 07/07/2015



Page: 1 of 2

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract. An electronic certificate can be authenticated online. Printed copies can be obtained at www.bsi.com/online/Default.htm

Information and Contact: BSI, Kitemark Court, Davy Avenue, King's Lynn, Norfolk, PE30 3BA, UK. Tel: +44 (0)1509 545454
BSI Registered UK Limited, registered in England under number 3855311 at 389 Chiswick High Road, London W4 4AL, UK.
A Member of the BSI Group of Companies.



bsi.

Certificate of Registration

QUALITY MANAGEMENT SYSTEM - ISO 13485:2003 & EN ISO 13485:2012

This is to certify that

Conmet, LLC
24/1 Onezhskaya str.
Moscow
125413
Russian Fed.

Holds Certificate Number: **MD 519196**

and operates a Quality Management System which complies with the requirements of ISO 13485:2003 & EN ISO 13485:2012 for the following scope:

The design and manufacture of Titanium implants and distraction devices for cranio-maxillo-facial surgery, Titanium Dental implant systems and Titanium spine systems.

For and on behalf of BSI:

Gary Fenton, Global Assurance Director

Originally registered: 15/08/2007

Latest Issue: 06/08/2013

Expiry Date: 24/08/2016



Page: 1 of 1

...making excellence a habit™

This certificate was issued electronically and remains the property of BSI and is bound by the conditions of contract. An electronic certificate can be authenticated online. Printed copies can be obtained at www.bsi.com/online/Default.htm

Information and Contact: BSI, Kitemark Court, Davy Avenue, King's Lynn, Norfolk, PE30 3BA, UK. Tel: +44 (0)1509 545454
BSI Registered UK Limited, registered in England under number 3855311 at 389 Chiswick High Road, London W4 4AL, UK.
A Member of the BSI Group of Companies.

bsi.

Certificate

Full Quality Assurance

No. **CE 519183**

Issued to:

Conmet, LLC
24/1 Onezhskaya str.
Moscow
125413
Russian Fed.

In respect of:

The design and manufacture of Titanium implants and distraction devices for cranio-maxillo-facial surgery, sterile and non-sterile Titanium Dental implant systems and Titanium spine systems.

on the basis of our examination under the requirements of Council Directive 93/42/EEC, Annex II, Section 3.2.

For and on behalf of BSI, a Notified Body for the above Directive (Notified Body Number 0066):

Gary Fenton, Global Assurance Director

First Issued: 3 Mar 2010

Date: 28 Aug 2013

Expiration Date: 3 Mar 2015

Page: 1 of 1


Conditions of Approval
Validity of this certificate is conditional on the quality system being maintained to the requirements of the Directive. This approval excludes all products designed and/or manufactured by a third party on behalf of the company named on this certificate, unless specifically agreed with BSI.

...making excellence a habit™


Information and Contact: BSI, Kitemark Court, Davy Avenue, King's Lynn, Norfolk, PE30 3BA, UK. Tel: +44 (0)1509 545454
BSI Registered UK Limited, registered in England under number 3855311 at 389 Chiswick High Road, London W4 4AL, UK. A member of the BSI Group of Companies.

1.5 Пластины сетчатые для краниопластики

толщина 0.6 мм




Пластины сетчатые плоские (100x100x0.6 мм)




Пластины сетчатые сферические (0.6 мм)

Радиусы сферы - R100, 130, 170 мм



R, мм	H, мм
100	24
130	18.5
170	14

1.5 Структура сетки



Ажур

517.03

Пластина сетчатая, плоская, ажур

517.13

Пластина сетчатая, сферическая, ажур, R100 мм

517.23

Пластина сетчатая, сферическая, ажур, R130 мм

517.43

Пластина сетчатая, сферическая, ажур, R170 мм


517.04

Пластина сетчатая, плоская, ажур, толщина 0.3 мм

1.5 Пластина сетчатая динамическая

толщина 0.6 мм


517.05




Пластина сетчатая, плоская, динамическая
(120x120x0.6 мм)

1.5 Винты самонарезающие для краниопластики

Кат. N	Диаметр, мм	Длина, мм
501.00	1.5	4
501.01	1.5	5
501.02	1.5	6




Кат. N	Диаметр, мм	Длина, мм
запасные		
502.00	1.8	5



1.5 Винты самосверлящие


Кат. N	Диаметр, мм	Длина, мм
501.10	1.5	4
501.11	1.5	5
501.12	1.5	6
501.13	1.5	8
501.14	1.5	10




1.5 Мини-пластины для дна глазницы

толщина 0.6 мм


M 1:1



15.04




15.02



15.03


1.5 Инструменты для краниопластики

522.00



Рукоятка поворотная универсальная
(Титан. сплав)

522.11




Жало сменное
(Нерж. сталь, лезвие с покрытием TiN)

519.05



Сверло перовое 1.2 x 5 мм
(Нерж. сталь)

521.01



Отвертка под угловой наконечник
(Нерж. сталь, лезвие с покрытием TiN)

2.0 Мини-пластины

толщина 1.0 мм

М 1:1



Прямая, длинная,
4 отв.

05.00



Прямая,
4 отв.

04.00



Прямая,
6 отв.

04.01



Прямая,
8 отв.

04.02




L-образная,
короткая, левая

06.00




L-образная,
короткая, правая

06.01



L-образная,
длинная, левая

07.00



L-образная,
длинная, правая

07.01