



HRCc[®] Cervical Locking Cage



HRCc® Cervical Locking Cage

- HRCc® is a Cervical Locking cage, supporting primary and secondary stability of the cage (thanks to its locking mechanism). The cage is inserted according to standard technique and then the surgeon locks the cage by simple rotation of the instrument.
- The blade penetrates into vertebral plates and prevents primary and/or secondary migration of the cage.
- May be used stand alone

Key Benefits:

- Easy surgical technique.
- Primary and secondary stability supported.
- X-Ray witness: Titanium blade and spot to verify the correct position of the cage.
- Stops on the instrument help to avoid accidental nerve compression.
- Excellent distribution of impacting shocks.
- All of the HRCc® Cervical Locking Cage are in PEEK Optima LT™ for a better biocompatibility and radiolucency.
- Easy Inventory.
- Large fusion surface.
- Excellent distribution of stress thanks to a large bearing surface to reduce risk of subsidence.
- Anatomical shape for a physiological fit. Sizes in 4, 5, 6, 7, 8*, 9*, normal and large*.
- Laser marking on the blade.
- In case of revision, the surgeon unlocks the cage by simple rotation of the instrument.
- Custom fit bone substitutes are available for each size (supplied separately).

Implants belong to class IIb
Instruments belong to class I
Trial cages belong to class IIa
Bone substitutes belong to class III
X-Rays : Journal Rachis, Vol. 17, n°1,
March 2005
LPP Code : for cage 8144129, for
bone substitute 8128024
CAT-HRCC-E3
Revision date : 27/04/2021



EUROSPINE
Life matters, We care.

Our vertebral column rests on pillars of innovation, the heart of our engagement beats stronger and stronger thanks to the experience accumulated since 1998.

Our spirit is marked by the concern for a unique meticulousness and Continuous simplicity in vertebral implant design.

We offer constantly new solutions that unravel the complexity of the spine, eases the fears of the patient and reinforce the expertise of the surgeon

Eurospine shares its vision around the entire world to ease the pain of the patient.

Our contribution makes the difference across the world because the unexpected could happen to anyone, and at any moment

Life Matters, We Care

Manufacturer : EUROSPINE Ltd
3 Avenue des Dahlias,
94240 L'Hay Les Roses, France
Tel : +33 (0)1 46 86 60 07
Fax : +33 (0)1 46 86 66 52
www.eurospine.com

*Optional (**) Up to the choice of surgeon
These images are none contractual