

ISIS™

PLIF PEEK Cage


Anatomic Design

Rigid fixation



Lordosis Maintenance

0, 4, 8 angles



Anti-backout

Surface design



Wide Window

Good bone ingrowth



Titanium Markers

Easy identification and positioning




Ordering Information

(Unit : mm)

| Article no. | Height | Width | Length | Angle | Article no. | Height | Width | Length | Angle |
|-------------|--------|-------|--------|-------|-------------|--------|-------|--------|-------|
| NP5080 | 8 | 11 | 24 | 0 | NP5114 | 11 | 11 | 24 | 4 |
| NP5090 | 9 | 11 | 24 | 0 | NP5124 | 12 | 11 | 24 | 4 |
| NP5100 | 10 | 11 | 24 | 0 | NP5134 | 13 | 11 | 24 | 4 |
| NP5110 | 11 | 11 | 24 | 0 | NP5144 | 14 | 11 | 24 | 4 |
| NP5120 | 12 | 11 | 24 | 0 | NP5154 | 15 | 11 | 24 | 4 |
| NP5130 | 13 | 11 | 24 | 0 | NP5108 | 10 | 11 | 24 | 8 |
| NP5084 | 8 | 11 | 24 | 4 | NP5118 | 11 | 11 | 24 | 8 |
| NP5094 | 9 | 11 | 24 | 4 | NP5128 | 12 | 11 | 24 | 8 |
| NP5104 | 10 | 11 | 24 | 4 | NP5138 | 13 | 11 | 24 | 8 |

ISIS™ PLIF PEEK Cage

User-friendly designed Instrument



| INSTRUMENT | |
|------------|--------------|
| ARTICLE NO | NAME |
| SCPI-001 | RETRACTOR 6 |
| SCPI-002 | RETRACTOR 11 |
| SCPI-020 | SHAVER 7 |
| SCPI-021 | SHAVER 8 |
| SCPI-022 | SHAVER 9 |
| SCPI-023 | SHAVER 10 |
| SCPI-024 | SHAVER 11 |
| SCPI-025 | SHAVER 12 |
| SCPI-026 | SHAVER 13 |
| SCPI-027 | SHAVER 14 |
| SCPI-031 | DISTRACTOR 8 |
| SCPI-032 | DISTRACTOR 9 |

| INSTRUMENT | |
|------------|--------------------|
| ARTICLE NO | NAME |
| SCPI-033 | DISTRACTOR 10 |
| SCPI-034 | DISTRACTOR 11 |
| SCPI-035 | DISTRACTOR 12 |
| SCPI-036 | DISTRACTOR 13 |
| SCPI-037 | DISTRACTOR 14 |
| SCPI-043 | HOLDER |
| SCPI-044 | SHEATH |
| SCPI-045 | IMPACTOR |
| SCPI-046 | BONE PACKING BLOCK |
| SCPI-047 | BONE PACKING BAR |
| SCPI-048 | CURETTE SQUARE |
| SCPI-051 | CONTAINER |
| SCPI-052 | IMPLANT CASE |