


INFORMATIONS CARD - A2L®

BLUE = REF OF IMPLANT GREEN= REF OF STERILE IMPLANTS ORANGE = REF OF INSTRUMENTS/TRIALS

REF	STERILE REFERENCES	Tr (mm)	t (mm)	L (mm)	P (mm)	CORRESPONDING TRIAL	LORDOSIS
A2L8-4	A2L8-4S	8	6	35.5	26	AIF8-4	4°
A2L9-4	A2L9-4S	9	7	35.5	26	AIF9-4	4°
A2L10-4	A2L10-4S	10	8	35.5	26	AIF10-4	4°
A2L11-4	A2L11-4S	11	9	35.5	26	AIF11-4	4°
A2L12-4	A2L12-4S	12	10	35.5	26	AIF12-4	4°
A2L13-4	A2L13-4S	13	11	35.5	26	AIF13-4	4°
A2L15-4	A2L15-4S	15	13	35.5	26	AIF15-4	4°
A2L9-9	A2L9-9S	9	5	35.5	26	AIF9-9	9°
A2L10-9	A2L10-9S	10	6	35.5	26	AIF10-9	9°
A2L11-9	A2L11-9S	11	7	35.5	26	AIF11-9	9°
A2L12-9	A2L12-9S	12	8	35.5	26	AIF12-9	9°
A2L13-9	A2L13-9S	13	9	35.5	26	AIF13-9	9°
A2L15-9	A2L15-9S	15	11	35.5	26	AIF15-9	9°
A2L10-13	A2L10-13S	10	4	35.5	26	AIF10-13	13°
A2L11-13	A2L11-13S	11	5	35.5	26	AIF11-13	13°
A2L12-13	A2L12-13S	12	6	35.5	26	AIF12-13	13°
A2L13-13	A2L13-13S	13	7	35.5	26	AIF13-13	13°
A2L15-13	A2L15-13S	15	9	35.5	26	AIF15-13	13°




AIP
HOLDER

AIR6
ROTATOR

INFORMATIONS CARD - A2L® Large

BLUE = IMPLANT REF GREEN = REF OF STERILE IMPLANTS ORANGE = REF OF INSTRUMENTS/TRIALS

	REF	STERILE REFERENCES	Tr (mm)	t (mm)	L (mm)	P (mm)	CORRESPONDING TRIAL	LORDOSIS
	A2LL8-4	A2LL8-4S	8	6	40	26	AIFL8-4	4°
	A2LL9-4	A2LL9-4S	9	7	40	26	AIFL9-4	4°
	A2LL10-4	A2LL10-4S	10	8	40	26	AIFL10-4	4°
	A2LL11-4	A2LL11-4S	11	9	40	26	AIFL11-4	4°
	A2LL12-4	A2LL12-4S	12	10	40	26	AIFL12-4	4°
	A2LL13-4	A2LL13-4S	13	11	40	26	AIFL13-4	4°
	A2LL15-4	A2LL15-4S	15	13	40	26	AIFL15-4	4°
	A2LL9-9	A2LL9-9S	9	5	40	26	AIFL9-9	9°
	A2LL10-9	A2LL10-9S	10	6	40	26	AIFL10-9	9°
	A2LL11-9	A2LL11-9S	11	7	40	26	AIFL11-9	9°
	A2LL12-9	A2LL12-9S	12	8	40	26	AIFL12-9	9°
	A2LL13-9	A2LL13-9S	13	9	40	26	AIFL13-9	9°
	A2LL15-9	A2LL15-9S	15	11	40	26	AIFL15-9	9°
	A2LL10-13	A2LL10-13S	10	4	40	26	AIFL10-13	13°
	A2LL11-13	A2LL11-13S	11	5	40	26	AIFL11-13	13°
A2LL12-13	A2LL12-13S	12	6	40	26	AIFL12-13	13°	
A2LL13-13	A2LL13-13S	13	7	40	26	AIFL13-13	13°	
A2LL15-13	A2LL15-13S	15	9	40	26	AIFL15-13	13°	
							AIP HOLDER AIR6 ROTATOR	