




EUROSPINE
Life matters, We care.

INFORMATIONS CARD - HRCC®

BLUE = REF. IMPLANT GREEN= REF. STERILE IMPLANYS ORANGE = REF. INSTRUMENTS/TRIALS

	REF	REF. STERILE	Tr (mm)	T (mm)	L (mm)	K (mm)	P (mm)	Corresponding Trial	LORDOSIS
	HRCC4SD	HRC4SD	4.2	3	16.6	12.8	14.4	HRCCIFS4	4.8°
	HRCC5SD	HRC5SD	5	4	16.6	12.8	14	HRCCIFS5	4°
	HRCC6SD	HRC6SD	6	5	16.6	12.8	14	HRCCIFS6	4°
	HRCC7SD	HRC7SD	7	6	16.6	12.8	14	HRCCIFS7	4°
	HRCC8SD	HRC8SD	8	7	16.6	12.8	14	HRCCIFS8	4°
	HRCC9SD	HRC9SD	9	8	16.6	12.8	14	HRCCIFS9	4°
	HRCC4LD	HRC4LD	4.2	3	18.6	14.8	14.4	HRCCIFL4	4.8°
	HRCC5LD	HRC5LD	5	4	18.6	14.8	14	HRCCIFL5	4°
HRCC6LD	HRC6LD	6	5	18.6	14.8	14	HRCCIFL6	4°	
HRCC7LD	HRC7LD	7	6	18.6	14.8	14	HRCCIFL7	4°	
HRCC8LD	HRC8LD	8	7	18.6	14.8	14	HRCCIFL8	4°	
HRCC9LD	HRC9LD	9	8	18.6	14.8	14	HRCCIFL9	4°	

	<p>HRCCIRS Rotateur</p> <p>HRCCIP Préhenseur</p>
--	--



EUROSPINE
Innovative Spinal Implants

CAGE INTERSOMATIQUE HRC
Lombaire Alif A2L Locking cage

CONDITIONS D'UTILISATION