

PediGuard[®] *Curv*



Product Presentation

The PediGuard[®] technology is now available in a curved tip.

Like the original PediGuard[®], the new PediGuard[®] Curv takes the place of the standard gearshift and facilitates pilot hole preparation for **accurate screw placement.**

Now surgeons have a choice.



PediGuard[®] Curv


PediGuard[®]

The new PediGuard[®] Curv offers the following benefits:

- **A tapered tip** to ease penetration through the pedicle *and* removal from the pedicle when drilling is complete.
- **A sense of directionality.** As surgeons become comfortable with the direction of the curve, we anticipate enhanced awareness of the location of the possible breach for proper redirection.
- Surgeons who prefer a curved gearshift will appreciate the **tactile feel** of the new PediGuard[®] Curv.




PediGuard[®] Curv




- Reference: **PI-AU450**
- Designation: **PediGuard[®] Curv**

Graduated diameter
from the tip to 30mm



Total length of 55mm

PediGuard[®] Curv



**Placing Navigation Back in Your Hands.
Now you have a choice.**