

Komponenten

Components

Implantate Implants	Länge x Breite Length x Width	Höhe Height	Winkel Angle	Artikelnummer Product number
pezo-P	24 x 10 mm	7-14 mm*	5°	CS 3301-07 to -14
pezo-P	29 x 10 mm	7-14 mm*	5°	CS 3302-07 to -14
pezo-P	29 x 10 mm	8-14 mm*	12°	CS 3303-08 to -14
pezo-T	29 x 11 mm	7-14 mm*	5°	CS 3315-07 to -14
pezo-T	29 x 14 mm	7-14 mm*	5°	CS 3316-07 to -14
pezo-T	34 x 16 mm	7-14 mm*	5°	CS 3317-07 to -14
pezo-A	28 x 22 mm	8-16 mm**	8°	CS 3331-08 to -16
pezo-A	31 x 25 mm	8-16 mm**	0°	CS 3332-08 to -16
pezo-A	31 x 25 mm	8-16 mm**	8°	CS 3333-08 to -16
pezo-A	31 x 25 mm	8-16 mm**	12°	CS 3334-08 to -16
pezo-A	34 x 28 mm	8-16 mm**	0°	CS 3335-08 to -16
pezo-A	34 x 28 mm	8-16 mm**	8°	CS 3336-08 to -16
pezo-A	34 x 28 mm	10-16 mm**	12°	CS 3337-10 to -16

* Standardbestückung für pezo-P und pezo-T bis Höhe 13 mm, Höhe 14 mm nur auf Anfrage

** 2 mm Schritte

* Standard equipment for pezo-P and pezo-T up to height 13 mm; height 14 mm upon request

** 2 mm increments